
## OBM Positioning Strips

### Instructions for Use:

Natus Medical Incorporated  
DBA Excel-Tech Ltd. (XLTEK)  
2568 Bristol Circle  
Oakville, Ontario, L6H 5S1 Canada  
Tel: +1 905-829-5300  
Website: [natus.com](http://natus.com)

EU Representative / Importer  
Natus Manufacturing Limited  
IDA Business Park  
Gort, Co. Galway, Ireland

Natus Nicolet UK Ltd  
Baynards Green Trading Estate  
Prospect House  
Oxfordshire  
Bicester  
OX27 7SG  
England, United Kingdom

Arazy Group Switzerland GmbH  
Bruderholzallee 53  
4059 Basel  
Switzerland  
[swiss.ar@arazygroup.com](mailto:swiss.ar@arazygroup.com)

Rx Only

### Associated product part numbers:

OBM00047

## English

### Description:

The OBM Positioning Strip is a single-use tool that uses anatomical landmarks on a neonatal head to mark positions for placement of amplitude-integrated EEG (aEEG) leads. Each package contains ten each of Term and Pre-Term sizes.

### Intended Use:

The intended use of the Olympic Brainz Monitor (OBM) positioning strips is to assist with the alignment of disposable surface recording electrodes.

### Intended User and Target Patient Group:

The product may be used with neonatal patients or infants (defined as birth to 28 days post-delivery, and corresponding to a post-conceptual age of 24 to 46 weeks), under the direction of a physician.

### Clinical Benefits:

Facilitates obtaining the recording of an EEG/PSG study to detect any irregularities indicative of various brain or sleep disorders.

### Contraindications and Side Effects:

There are no known contraindications or side effects for marking performed with the OBM Positioning Strip. See Warnings or Precautions below.

### Operating Instructions:

- Follow pictorial and written instructions on the front and back of the OBM positioning strip.
- Mark indicated locations for placement of aEEG leads.
- Dispose of the paper strip following a single instance of use.

### Understanding Warnings and Cautions Statements:

 <b>CAUTION</b>
<b>Refers to a hazardous situation that could result in minor or moderate injury or material damage if not avoided.</b>
<ul style="list-style-type: none"> <li>Information on how the hazardous situation is avoided.</li> </ul>

### Warnings and Precautions:

 <b>CAUTION</b>
<b>Device dropped or damaged in transit/use could lead to loss of function or delayed diagnosis.</b>
<ul style="list-style-type: none"> <li>Inspect the device prior to each use and do not use if damaged.</li> </ul>

 <b>CAUTION</b>
<b>Device reused on another patient leads to cross infection or loss of performance.</b>
<ul style="list-style-type: none"> <li>Do not reuse or sterilize.</li> </ul>
<b>Improper usage of OBM positioning strip might lead to wrong placement.</b>
<ul style="list-style-type: none"> <li>Only clinical personnel with clinical expertise in aEEG and EEG lead placement should utilize the OBM Positioning Strip.</li> </ul>

### Environmental Specifications:

#### Operating Conditions:

- Temperature: +10°C (+50°F) to +30°C (+86°F)
- Relative Humidity: 20% to 80%
- Pressure: 70 kPa to 106 kPa

#### Storage Conditions:

- Temperature: 0°C (+32°F) to +30°C (+86°F)
- Relative Humidity: 20% to 80%
- Pressure: 50 kPa to 106 kPa

### Compliance Standards:

- ISO 10993-1: 2018 Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process
- ETS 300 019-2-1 Environmental Engineering (EE); Environmental conditions and environmental tests for telecommunications equipment; Part 2-1: Specification of environmental tests; Storage
- ETS 300 019-2-2 Environmental Engineering (EE); Environmental conditions and environmental tests for telecommunications equipment; Part 2-2: Specification of environmental tests; Transportation
- ASTM D4169-16 Standard Practice for Performance Testing of Shipping Containers and Systems for Vibration


### Disposal Instructions:

Following use, dispose of adhesive electrodes with infectious waste.

### Disclaimer:

Natus Medical Incorporated DBA Excel-Tech Ltd. (Xltek) is not responsible for injury, infection or other damage resulting from the use of this product. Any serious incident that has occurred in relation to the device should be reported to Natus Medical Incorporated DBA Excel-Tech Ltd. (Xltek) and the competent authority of the Member State in which the user and/or patient is established. Refer to the Natus website for an electronic copy of this document.

## Glossary of Symbols

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title as per Reference Standard	Explanation
	-	-	An indication of Medical device	This product is a medical device.
	21 CFR Part 801.109(b)(1)	Labeling-Prescription devices.	Prescription only	Indicates the product is authorized for sale by or on the order of a licensed healthcare practitioner.
	ISO 15223-1 Symbol 5.1.1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Manufacturer	Indicates the medical device manufacturer.
	ISO 15223-1 Symbol 5.1.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Authorized representative in the European Community	Indicates the Authorized representative in the European Community.
	ISO 15223-1 Symbol 5.1.6	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1 Symbol 5.1.5	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Batch or Lot code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1 Symbol 5.1.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Date of manufacture	Indicates the date when the medical device was manufactured.
	ISO 15223-1 Symbol 5.4.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	ISO 60601-1 Table D.1 #11	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.	Operating instructions	

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title as per Reference Standard	Explanation
	ISO 15223-1 Symbol 5.4.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Caution: Read all warnings and precautions in instructions for use	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	ISO 60601-1 Table D.1 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.		
	ISO 60601-1 Table D.2 #2	Medical electrical equipment — Part 1: General requirements for basic safety and Essential performance.	General warning sign	Indicates a hazard of potential personal injury to patient or operator.
	ISO 15223-1 Symbol 5.2.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1 Symbol 5.3.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Temperature limit	Indicates the (storage) temperature limits to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.3.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Humidity limitation	Indicates the range of (storage) humidity to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.4.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	ISO 15223-1 Symbol 5.4.5 (Reference Annex B for the general prohibition symbol)	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Not made with Natural Rubber Latex	Indicates a medical device that is not made with natural rubber latex.
	2012/19/EU	Waste Electrical and Electronic Equipment (WEEE)	Disposal at end of operating life instructions	Indicates that electrical and electronic equipment waste should not be discarded together with unseparated waste but must be collected separately.
	ISO 15223-1 Symbol 5.3.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Keep away from sunlight	Indicates a medical device that needs protection.
	Swiss Medical Device Ordinance (MedDO)	Swiss Medical Device (MedDO).	Indicates the Authorized Representative in Switzerland	Indicates the Authorized Representative in Switzerland.
	UKCA Medical Device Regulation (SI 2022 No 618, as amended) (UK MDR 2002)	UKCA Medical Device Regulation.	UKCA Mark	In compliance with the United Kingdom technical conformity.

Product Part Numbers	Product Description
OBM00047	Pre-Term Size Strips, 10 each

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