# natus®

## Disposable Snap Electrodes

Instructions for Use

Brugsanvisning

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Mode d'emploi

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## English

### Description:

Snap Electrodes are used in conjunction with Snap Leads and amplifiers for monitoring electromyogenic and physiologic data signals commonly recorded from chest, extremities and faces of patients.

### Intended Use:

Evoked potential Disposable Snap Electrodes are intended to be used with compatible 1.5 mm touch proof snap electrode leads for testing of the electrical activity from the brain, spinal nerves, or sensory receptors in response to specific external stimulation.

### Intended User and Target Patient Group:

The product is for use only under the direction of trained medical professionals in acquiring electrophysiologic signals. For use on Infants.

### Clinical Benefits:

Facilitates obtaining the recording of an EEG/PSG study to detect any irregularities indicative of various brain or sleep disorders.

### Contraindications and Side Effects:

There are no known contraindications or side effects for procedures performed with Disposable Snap Electrodes. See Warnings or Precautions below.

### Operating Instructions:

- Open the pouch just before use.
- Prepare the electrode site prior to application.
- Connect the appropriate snap lead to the snap connector on the electrode.
- Remove the electrode from its backing.
- Apply to the patient making sure to place even pressure across the entire electrode area.
- Electrodes can be used to record or stimulate.
- Re-seal pouch after each use.
- Once recording is complete, detach and discard the disposable snap electrode.

### Understanding Cautions Statements:

<b>⚠ CAUTION</b>
<b>Refers to a hazardous situation that could result in minor or moderate injury or material damage if not avoided.</b>
<ul style="list-style-type: none"><li>• Information on how the hazardous situation is avoided.</li></ul>

### Warnings and Precautions:

<b>⚠ CAUTION</b>
<b>Device dropped or damaged in transit/use could lead to loss of function or delayed diagnosis.</b>
<ul style="list-style-type: none"><li>• Inspect the device prior to each use and do not use if damaged.</li></ul>
<b>Electrode reused on another patient leads to cross infection or loss of performance.</b>
<ul style="list-style-type: none"><li>• Do not reuse or sterilize.</li></ul>
<b>Data cannot be recorded due to improper placement leading to delay in procedure.</b>
<ul style="list-style-type: none"><li>• Assure proper connections and check signal quality before using.</li></ul>
<b>Electrodes left in place for prolonged periods of time may be associated with skin redness.</b>
<ul style="list-style-type: none"><li>• Do not use the electrodes for prolonged periods of time.</li></ul>

### Repositioning of electrode can degrade the signal quality.

- Electrode should be replaced if the adhesive no longer provides a uniform connection to the skin surface.

### Modifications to the device can affect function and performance.

- Do not modify this equipment without authorization of the manufacturer.

### Environmental Specifications:

#### Operating Conditions:

- Temperature: +10°C (+50°F) to +30°C (+86°F)
- Relative Humidity: 20% to 80%
- Pressure: 70 kPa to 106 kPa

#### Storage Conditions:

- Temperature: 0°C (+32°F) to +30°C (+86°F)
- Relative Humidity: 20% to 80%
- Pressure: 50 kPa to 106 kPa

### Compliance Standards:

- ISO 10993-1:2018 Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process
- ETS 300 019-2-1 Environmental Engineering (EE); Environmental conditions and environmental tests for telecommunication equipment; Part 2-1: Specification of environmental tests; Storage
- ETS 300 019-2-2 Environmental Engineering (EE); Environmental conditions and environmental tests for telecommunications equipment; Part 2-2: Specification of environmental tests; Transportation
- ASTM D4169-16 Standard Practice for Performance Testing of Shipping Containers and Systems for Vibration
- EN 60601-1:2006/A1:2013 and IEC 60601-1:2005/A1:2012 Ed. 3.1 clause 8.5.2.3

### Disposal Instructions:

Following use, dispose of adhesive electrodes with infectious waste.

### Disclaimer:

Natus Medical Denmark Aps is not responsible for injury, infection or other damage resulting from the use of this product.

Any serious incident that has occurred in relation to the device should be reported to Natus Manufacturing Limited and the competent authority of the Member State in which the user and/or patient is established.

Refer to the Natus website for an electronic copy of this document.

## Glossary of Symbols

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title as per Reference Standard	Explanation
<b>Medical Device</b>	-	-	An indication of Medical device	This product is a medical device.
<b>Rx Only</b>	21 CFR Part 801.109(b)(1)	Labeling-Prescription devices.	Prescription only	Indicates the product is authorized for sale by or on the order of a licensed healthcare practitioner.
	ISO 15223-1 Symbol 5.4.5 (Reference Annex B for the general prohibition symbol)	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Not made with Natural Rubber Latex	Indicates a medical device that is not made with natural rubber latex.
	-	-	Quantity	Number of parts in a package.
	ISO 15223-1 Symbol 5.1.1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Manufacturer	Indicates the medical device manufacturer.
	ISO 15223-1 Symbol 5.1.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Date of manufacture	Indicates the date when the medical device was manufactured.
	ISO 15223-1 Symbol 5.1.5	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Batch or Lot code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1 Symbol 5.1.6	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1 Symbol 5.4.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	ISO 60601-1 Table D.1 #11	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.	Operating instructions	
	ISO 15223-1 Symbol 5.4.3 Annex A #A.15	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Consult instructions for use	Indicates an instruction to consult an electronic instructions for use (eIFU).
	ISO 60601-1 Table D.2 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.	Follow instructions for use	Refer to instruction manual/ Booklet.  NOTE on ME EQUIPMENT "Follow instructions for use"
	ISO 15223-1 Symbol 5.4.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Caution: Read all warnings and precautions in instructions for use	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	ISO 60601-1 Table D.1 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.		
	ISO 15223-1 Symbol 5.1.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Use-by date	Indicates the date after which the medical device is not to be used.  Note: This symbol shall be accompanied by a date to indicate that the medical device should not be used after the end of the year, month or day shown.
	ISO 15223-1 Symbol 5.2.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1 Symbol 5.3.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Temperature limit	Indicates the (storage) temperature limits to which the medical device can be safely exposed.

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title as per Reference Standard	Explanation
	ISO 15223-1 Symbol 5.3.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Humidity limitation	Indicates the range of (storage) humidity to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.4.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.

Product Part Numbers	Product Description
8-64-21000	Disposable Snap Electrode 3/Pkg.
8-64-21800	Disposable Snap Electrode 25 x 3/pkg (75 total)
8-64-21600	Disposable Wet Gel Snap Electrode 20/pkg
8-64-21601	Disposable Adhesive Electrode 10/pkg

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