


DURAFORM[®] Dural Graft Implant

Designed for Superior Performance

The DURAFORM Dural Graft Implant is the solution for the difficult handling of dural substitutes. The unique design and manufacturing process provide an implant with better handling characteristics and greater tear resistance than other onlay dural graft materials. The DURAFORM provides superior leak resistance with two different textured surfaces for enhanced wet handling capabilities and excellent anatomical conformity.


Features:

- Textured top surface makes the DURAFORM implant easy to grip and manipulate, by hand or with forceps
- Soft, smooth bottom side that gently molds to the complex surface structures of the brain and dura
- Offers excellent anatomical conformity to help prevent CSF leakage
- Collagen-based biocompatible material exhibits superior tensile strength and wet handling capabilities
- The strength of the DURAFORM Dural Graft allows for improved manipulation of the dural substitute
- Extended shelf life (2 years)



Enhanced Confidence with DURAFORM

- CSF leak resistance is a key factor in the overall effectiveness and value of dural substitutes
- Made of a collagen-based biocompatible material obtained from a geographical BSE risk level 1 source
- When wet, the DURAFORM implant is stronger and handles better than other onlay graft dural substitutes.*



Ordering Information

Item Number	Size	Units/Case
801480US	4 in x 5 in (10 cm x 12.5 cm)	1
801478US	3 in x 3 in (7.5 cm x 7.5 cm)	5
801477US	3 in x 3 in (7.5 cm x 7.5 cm)	1
801476US	2 in x 2 in (5 cm x 5 cm)	5
801475US	2 in x 2 in (5 cm x 5 cm)	1
801474US	1 in x 3 in (2.5 cm x 7.5 cm)	5
801473US	1 in x 3 in (2.5 cm x 7.5 cm)	1
801472US	1 in x 1 in (2.5 cm x 2.5 cm)	5

Indications: The DURAFORM Dural Graft Implant is intended for use in procedures where the repair or substitution of the patient's dura mater is needed.

Contraindications: The DURAFORM Dural Graft Implant is contraindicated in patients with known sensitivity to bovine-derived materials.

Do not use to repair spinal neural tube defects or for anterior spinal surgery with dural resection.

Use with caution in infected areas.

Natus does not recommend use to cover large defects at the skull base following surgery or dural defects involving mastoid air cells.



Natus Neuro
 3150 Pleasant View Road
 Middleton, WI 53562 USA
 Tel: 1-800-356-0007
 1-608-829-8500
 Fax: 1-608-829-8709
www.natus.com

To learn more about Natus Neuro products

Call: **1-800-356-0007**

Visit: **www.natus.com**

*Data on file.

©2017 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.