

## Nicolet Cortical Stimulator for Neurosurgery


### Functional brain mapping for brain tumor resection

The Nicolet® Cortical Stimulator is a compact, digital and portable unit that aids in identification and mapping of the areas of the brain that control speech, sensory, vision and motor skills. The resulting map aids physicians and surgeons in planning surgical approach and procedure.

In addition to providing direct, constant current stimulation to the brain, the Nicolet Cortical Stimulator offers full control over pulse frequency and duration, train duration and current levels.

#### Features and Benefits:

- Displays selected and actual current delivered
- Ictal Disrupt feature for management of "after discharges"
- Small, compact size
- Easy to use controls
- Disposable probe
- Optional integrated software




# The new standard of care

The innovative ergonomic compact design, easy to control settings with disposable probe and enhanced functionality helps neurosurgeons when performing brain mapping procedures in the operating room during brain surgery, including tumor and epilepsy resections.

The Delivered Stim feature displays the actual stimulation delivered to the brain, allowing better control of stimulus levels. The stimulator can be used in the single pulse or continuous mode and shows a countdown of the remaining stimulation time as well as an active stimulation indicator light. Users have told us this device provides a new standard of care at their fingertips.

The Nicolet Cortical Stimulator received 510(K) clearance for the intended use of performing direct cortical stimulation as a stand alone device or as interfaced with its Nicolet LTM system through use of the related cortical stimulation software. The device is available both in the U.S. and internationally.



## Service

Excellence in customer and technical service is a core value of Natus Neurology.

### Essential Support Elements:

- Accessible and effective Technical Support
- Definitive technical documentation and knowledgeable installation teams
- Excellent replacement unit or spare part availability
- Extended warranty and service coverage programs
- Comprehensive, flexible customer training courses

## Supplies

Natus Neurology offers a full range of neurodiagnostic supplies, featuring:

- Dedicated and knowledgeable customer service
- Streamlined order processing to save you time
- Convenient online portal (US Customers only)  
[Natus Neuro Store – www.natusneurostore.com](http://www.natusneurostore.com)

To learn more about Natus Neurology Service Programs or our full line of Supplies, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**



### Natus Neurology Incorporated

3150 Pleasant View Road  
Middleton, WI 53562 USA  
Tel: 1-800-356-0007  
1-608-829-8500  
Fax: 1-608-829-8709  
[www.natus.com](http://www.natus.com)

©2016 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.