

December 10, 2018

Re: New Software Release

Dear Valued Customer:

Announcing the Release of NeuroWorks/SleepWorks 9.1 GMA 1

We are pleased to announce the release of a new software version for NeuroWorks and SleepWorks. NeuroWorks/SleepWorks version 9.1 GMA 1 is an enhancement that includes usability improvements containing but not limited to:

1. *Introduction of an LTM Study as a new NeuroWorks Study Type for EEG and Slideshow navigation tool*

An 'LTM Study' is a new study type introduced to act as a collection of NeuroWorks studies for a multi-day patient stay. In combination with Slideshow, it drastically simplifies review, reporting and pruning for multi-day studies.

2. *Mandatory Fields*

Allows configuring certain fields as mandatory during new or returning study operation.

3. *Improved Visual Montage Editor*

Visual Montage Editor adds the following new features:

- Creation of a standard Contralateral or Ipsilateral montage (based on 10-20 layout).
- Ability to create nearest neighbor montage (based on 10-20 layout).
- Ability to create a sequence of channels into a list of channel pairs.

4. Improved video playback

Added the ability to quickly page through the study with video by displaying a video image for every EEG page.

5. Reporting improvements

- **Automatically export study snapshots into reports**
Added the ability to setup report templates to automatically include EEG snapshots saved during study review
- **Export pre-defined note types to reports**
Added the ability to setup report templates to automatically include specific types of annotations

All customers who currently have NeuroWorks/SleepWorks version 9.0 are entitled to a free software update to 9.1 GMA 1. If you own one of the older NeuroWorks/SleepWorks versions you may order NeuroWorks/SleepWorks 9.1 GMA 1 as an upgrade.

For download links please contact our Technical support department at 1-800 387-7516 or email Oakville_Technical_Service@natus.com to make arrangements if you choose.

Thank you for being our healthcare partner!

Sincerely,

Gilad Shulkin
Director of EEG Marketing
Natus Neuro