



## Viking® on UltraPro S100

Continuing Excellence in  
Electrodiagnostic Systems

### Reliability, Performance, Portability – Exceptional Value

With over 60 years of proven experience, Natus Neuro introduces the Viking UltraPro S100. Featuring automatic reporting and user-friendly hardware and software, the Viking UltraPro S100 bundles high performance, ease of use, and exceptional value in a compact and reliable package.

#### Proven Reliability

Built on the 60-year heritage of Nicolet®, Dantec™, Xltek® and Teca® brands, Viking UltraPro S100 was developed in close cooperation with physicians and technologists and continues a tradition of providing reliable and dependable electrodiagnostic equipment.

#### High Quality

All-in-one electrodiagnostic system with dedicated controls, integrated stimulators and exceptional audio and signal quality.

#### Exceptional Value

Affordable equipment that is flexible and powerful enough to meet your diagnostic needs without compromising ease of use or quality.

#### Trouble-Free Installation

With only 3 cables, setup is easy and no installation support is required. Minimal set up time allows you to start performing studies immediately.

#### Ease of Use


The Viking UltraPro S100 allows you to efficiently complete procedures and examinations, providing industry leading quantitative data, analysis, and results quickly. The Viking UltraPro S100 also includes a robust Help menu to provide immediate support assistance for many areas, as well as industry leading instructional videos for any training follow-up or refresher.

#### Available with Evoked Potentials


Add any or all of the following EP's as an option for your Viking UltraPro S100.

SKU #	Description
828-061600	EP (3 Modality – SEP, VEP, AEP)
828-061500	Somatosensory EP (SEP)
828-061400	Visual EP (VEP)
828-061300	Auditory EP (AEP)


# Versatile Viking UltraPro S100 Software for Electrodiagnostic and Nerve Conduction Studies



Simple Setup



Electrodiagnostics with Scoring Tool




Automatic Summary Generation

## Simple workflow

The one-button EMR interface, distinctive function keys, online summaries and built-in reference value database saves time, generating reports faster, improving service to referring physicians and connecting more efficiently with electronic medical records.

## Focus on the patient, not the equipment

The UltraPro S100 is an all-in-one lightweight system that travels easily in a roller case. The available triple programmable footswitch supports hands-free operation. Three- and four-channel amplifiers allow electrodes to remain in place, accommodating the most common modalities (motor and sensory NCS, needle EMG and EP's), reducing test time and improving the patient experience. The UltraPro S100 is also available as a mobile cart-based system.



## Service & Support

A global service organization here for you

### Essential Support Elements:

- Accessible and effective Technical Support
- Definitive technical documentation and knowledgeable installation teams
- Excellent replacement unit or spare part availability
- Extended warranty and service coverage programs

## Supplies

Natus Neuro offers a full range of neurodiagnostic supplies.

- Dedicated and knowledgeable customer service
- Streamlined order processing to save you time
- Convenient online portal (US Customers only)  
Natus Neuro Store – [www.natusneurostore.com](http://www.natusneurostore.com)

To learn more about Natus Neuro Service Programs or our full line of Supplies, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**



Natus Neuro

3150 Pleasant View Road  
Middleton, WI 53562 USA  
Tel: 1-800-356-0007  
1-608-829-8500  
Fax: 1-608-829-8709

[neuro.natus.com](http://neuro.natus.com)