



## Pocket-Dop™ II

### Obstetric and Vascular Screening and Monitoring

#### Obstetric & Vascular Versatility

Pocket-Dop II is a full-featured Doppler with interchangeable obstetric and vascular probes. Use the 2 or 3 MHz obstetric probes to monitor the fetal heartbeat. Use the 5 or 8 MHz vascular probes to screen patients at risk for stroke, peripheral arterial disease and other occlusive vascular conditions or to monitor difficult pulses and obtain systolic blood pressures quickly and easily.

#### Obstetric Examination

The Pocket-Dop II is a valuable aid in detecting pregnancy and monitoring the fetal heartbeat. The built-in speaker amplifies fetal sounds to reassure the expectant mother. The 2 MHz waterproof obstetric probe is suited for underwater labor and delivery.

#### Vascular Applications

Use the Pocket-Dop II to compare simple Doppler blood pressures in the arm and ankle to help diagnose arterial disease in the lower extremities or quickly screen for lower extremity vascular conditions such as valvular incompetence.

#### Emergencies & Challenging Clinical Situations

The Pocket-Dop II can be used to accurately monitor pulses and systolic blood pressures, even in noisy environments; helps confirm blood flow to wound sites, locate vessels for catheterization, withdrawal, injection or IV therapy; clearly hear low systolic blood pressures with patients who are diabetic, elderly or in shock.



## Interchangeable Probes



### 2 MHz obstetric probe

Deep signal penetration and best used to determine fetal viability in larger women and later in pregnancy. It will consistently pick up the fetal heartbeat at 12 weeks and preferred for use during labor and delivery. The 2 MHz waterproof obstetric probe is suited for underwater labor and delivery.



### 3 MHz obstetric probe

Shallow penetration and more sensitive to early fetal heartbeats. It can pick up fetal signals as early as 8-10 weeks.



### 5 MHz vascular probe

Deeper penetration than the 8 MHz probe. It is most commonly used for post-surgical patients, locating hard-to-find pulses and monitoring deep veins and arteries in legs; and for ABI studies in larger patients.



### 8 MHz vascular probe

Shallow signal penetration and best used to assess superficial vessels. It is most commonly used for ankle/brachial index (ABI) studies, finding pedal pulses, wound care applications and locating hard-to-find pulses.

## System Configurations

Product	Catalog #
<b>Pocket-Dop II with</b>	
2 MHz obstetric probe .....	P200
3 MHz obstetric probe .....	P300
5 MHz vascular probe .....	P500
8 MHz vascular probe .....	P800
Two 2 MHz obstetric probes .....	P220
2 & 3 MHz obstetric probes .....	P230
2 MHz obstetric & 5 MHz vascular probes .....	P250
2 MHz obstetric & 8 MHz vascular probes .....	P280
Two 3 MHz obstetric probes .....	P330
3 MHz obstetric & 5 MHz vascular probes .....	P350
3 MHz obstetric & 8 MHz vascular probes .....	P380
Two 5 MHz vascular probes .....	P550
5 & 8 MHz vascular probes .....	P580
Two 8 MHz vascular probes .....	P880
2 MHz waterproof obstetric probe .....	P2WP
2 MHz waterproof and 3 MHz obstetric probes .....	P2W3
2 & 3 MHz obstetric and 5 & 8 MHz vascular probes .....	P2358
No probe (with accessories).....	PX00
No probe and no accessories .....	XPA002

## Technical Specifications

Weight .....	0.285 kg (10 oz.)
Dimension .....	3.17 x 5.71 x 10.79 cm (1.25 x 2.25 x 4.25 in)
Power Source .....	2 replaceable AA NiMH batteries
	110 volt or 220 volt battery charger

## Warranty against manufacture defects

Unit .....	10 years parts and 1 year labor against manufacture defects
Probes .....	1 year parts and labor against manufacture defects



### Natus Neuro

3150 Pleasant View Road  
Middleton, WI 53562 USA  
Tel: 1-800-356-0007  
1-608-829-8500  
Fax: 1-608-829-8709  
[neuro.natus.com](http://neuro.natus.com)

### Natus – Nicolet brand products

720 Corporate Circle, Suite A  
Golden, CO 80401 USA  
Tel: 1-877-842-7970  
1-303-431-9400  
Fax: 1-303-431-0429  
1-800-939-4639

## International System Configurations (220 volt for outside the US only)

Product	Catalog #
<b>Pocket-Dop II System with</b>	
2 MHz obstetric probe .....	P200I
3 MHz obstetric probe .....	P300I
5 MHz vascular probe .....	P500I
8 MHz vascular probe .....	P800I
2 and 3 MHz obstetric probes .....	P230I
2 MHz obstetric and 5 MHz vascular probes .....	P250I
2 MHz obstetric and 8 MHz vascular probes .....	P280I
3 MHz obstetric and 5 MHz vascular probes .....	P350I
3 MHz obstetric and 8 MHz vascular probes .....	P380I
Two 8 MHz vascular probes .....	P880I
5 and 8 MHz vascular probes .....	P580I
2 MHz waterproof obstetric probe .....	P2WPI
No probe and no accessories .....	XPA002

## Accessories

Product	Catalog #
<b>Probes</b>	
2 MHz obstetric probe .....	T200
3 MHz obstetric probe .....	T300
5 MHz vascular probe .....	T500
8 MHz vascular probe .....	T800
2 MHz waterproof probe .....	WP20
<b>Rollstands and Rollstand Accessories</b>	
With the bracket .....	ST1
With bracket & storage basket .....	ST2
Storage basket for rollstand .....	BSKT
Locking bracket only .....	B100
<b>Chargers</b>	
110 volt .....	C610
220 volt .....	C619I
Battery .....	C620

For more information, please contact  
your Nicolet Distributor or call:

**Toll Free: 877-842-7970**  
**International: 303-431-9400**  
**Fax: 303-431-0429**