

July 1, 2019

Dear Valued Customer,

Re: Protektor32 – EPWorks v6.1.2 software update

We are pleased to announce the release of the EPWorks v6.1.2 software update for Protektor32 IOM systems. Below are the high-level changes in version 6.1.2:

- Improved functionality within the Timeline to show status and ability to advance to any time position within the study.
- Improved the ability to advance EEG data when reviewing a study.
- Improved the configuration of the sound output for TIPS in BAEP monitoring.

Protektor32 customers with EPWorks v6.0 and v6.1 are entitled to the update to version 6.1.2 at no charge. To download this software free-of-charge, please contact our Technical Support department for further information. Please include your Protektor32 model and serial number in any correspondence.

Natus Technical Support Phone: 1-800-356-0007 or 1-608-829-8500

Natus Technical Support Email: OTS@natus.com

Natus Neuro is committed to developing the very best Neurodiagnostic products.

Thank you for being our healthcare partner and please feel free to contact your sales representative for further information.

Thank you,

Marketing
Natus Medical Incorporated
natus.