

August 2019

Dear Valued Customer,

**SUBJECT: Release of Nicolet EMG Software Version 22.3.0 Update for EDX, VikingQuest and UltraPro**

We are pleased to announce the release of the Nicolet EMG software v22.3.0 update for EDX, VikingQuest and UltraPro. Below are the high-level changes in v22.3.0:

- **Higher sampling rate at 48 kHz** for acquired and saved SNCS/MNCS data, providing high resolution data that improves data analysis precision.
- **Improved accuracy with AANEM Reference Values** with additional gender, BMI and height norms. Additionally, norms can be easily defined for pediatric populations.
- **Improved usability and workflow** with new software prompts improves ease of use, streamlines patient testing, and provides simple test sorting and easy merging/moving of tests. In-test benefits include the NCS marker amplitude/latency separation and convenient EP manual marking during data acquisition with real-time updates of the results table.
- **Improved dependability and performance in networked environments**, increasing speed and performance.

To obtain this software, please contact your Natus sales representative for further information. Please include your EDX, VikingQuest or UltraPro model and serial number in any correspondence.

**Software Orders:**

For software updates from V22.0 or later, order:

**Part# 688-628600 - Software Update to current version**

For software upgrades from V21.1 or earlier, order:

**Part# 688-628500 - Software Upgrade to current version**

For Support: Natus Technical Support Phone: 1-800-356-0007 or 1-608-829-8500  
Natus Technical Support Email: [MTS@natus.com](mailto:MTS@natus.com)

Natus is committed to developing the very best Neurodiagnostic products. Thank you for being our healthcare partner. We encourage you to contact your sales representative for further information on v22.3.0.

Thank you,

*Product Marketing*

---

Natus Medical Incorporated [www.natus.com](http://www.natus.com)

