

December 23rd 2019

Dear Valued Customer:

SUBJECT: InVisus® Ultrasound v6.2.8 software update

We are pleased to announce the release of the InVisus Ultrasound v6.2.8 software updates.

Here is a summary of the high-level improvements:

- **Bluetooth File Transfer Capability** – Makes it possible to wirelessly share images and video files to the paired computer system.
- **HDMI Output Option Selection** – Allows the user to select from HDMI display options of “Extend” and “Clone”. The new “Clone” option mirrors the InVisus display onto the external monitor.
- **Contrast Adjustment** – It is now possible to adjust Contrast (white/black ratio) while scanning.
- **New “Needle Preset”** – Highlights the needle during ultrasound-guided injection procedures.
- **Image Pan Deactivate Option** – Deactivating this display feature prevents the ultrasound image from being moved around within the screen. It remains in a fixed location.
- **Extended Temperature Feature** – This setting allows InVisus to scan longer during high temperatures, reducing the occurrence of overheating messages.
- **Pulse Wave (PW) PRF Adjustment** – Adjustments to the PRF are done by using the slide cursor.
- **Color Inversion and PRF for Color Doppler Mode** – Easily invert the colors by pressing and holding the Gain cursor button for 5 seconds. Adjusting the pulse repetition frequency (PRF) flow detection threshold allows the user to detect slower flows.
- **Manual measurements for the Pulsed Waved Doppler mode** – This function is used to manually measure the maximum and minimum speed of the flow.

InVisus customers are entitled to the update to version 6.2.8 at no charge. To download this software free-of-charge, please contact our [Technical Support](#) department for further information. Please include your InVisus model and serial number in any correspondence.

Natus Technical Support Phone: 1-800-356-0007 or 1-608-829-8500

Natus Technical Support Email: MTS@natus.com

Thank you for being a Natus healthcare partner, we are committed to developing the best Neurodiagnostic solutions in the industry. Please feel free to contact your sales representative for further information.

Thank you,

Marketing Team

Natus Medical Incorporated

Marketing@natus.com