



Reliable. Proven. Trusted.

The Natus Camino ICP Monitor (CAM02) is an advanced, touchscreen monitoring system that captures real-time data at the time of catheter insertion and trend data over a five-day period. Camino is a compact and portable device that provides tools for continuously determining and monitoring intracranial pressure (ICP) and intracranial temperature (ICT).

Advanced and versatile features

- Intuitive, easy-to-use menu system puts five days of patient information at your fingertips with rapid data export to USB
- Widescreen display ensures that critical data can be read from 10 feet away and up to a 30 degree offset from the monitor
- Ergonomic handle and convenient size for bedside pole attachment
- Tunneled or bolted catheter options



Camino Family of Catheters

Camino Flex (VTUN)

Therapeutic ventricular drainage of CSF while performing continuous ICP monitoring with the Camino Monitor

B-Series (1104B, 1104BT)

Parenchymal bolt catheters for measuring ICP and ICT

HM-Series (1104HM, 1104HMT)

Ventricular bolt catheters for measuring ICP and ICT

1104G

Subdural, post craniotomy catheter for measuring ICP


1104L

Parenchymal bolt catheter that secures through the Licox® bolted system for measuring ICP and O₂

Detailed information about the Natus Camino Catheters can be found in our product catalog.

Key Features of Camino catheters

- Easy and reliable way to measure ICP and ICT
- Camino bolt is quick and easy to place – Parenchymal and ventricular options available
- Accurate monitoring at the source
- Continuous parenchymal or ventricular monitoring



Versatility in draining and monitoring with a family of bolted and tunneled catheters

Natus Camino catheters measure their respective values at the tip of the catheter, providing continual ICP measurements. This design eliminates the need for a fluid filled system to communicate pressure to an external transducer.

Service

At Natus, we strive for excellence in customer and technical service.

Here's how we can help:

- Accessible and effective Technical Support
- Replacement unit and spare part availability
- Extended warranty and service coverage programs
- Comprehensive, flexible customer training courses

To learn more about Natus products, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. Licox is a registered trademark of Integra LifeSciences Corporation. P/N 023154D

natus

Natus Medical Incorporated

natus.com