


# Grass® Disposable Deep Cup Electrodes



## Trusted, Reliable, Convenient

### New 25-Packs!

#### Materials

##### Electrode

Silver/Silver Chloride coated ABS plastic

##### Cup Dimensions

Cup diameter ..... 10mm

Cup height ..... 3.0mm

Interior cup diameter ..... 6.7mm

##### Lead Wire

Core ..... Tin plated Copper

Jacket ..... PVC

##### Touchproof connector, 1.5mm female (DIN 42 802)

Connector ..... Tin plated Brass

Insulation ..... PVC

##### RoHS compliant

To learn more about Natus products, contact your local distributor or sales representative.

US customers call:  
**1-800-356-0007**

International customers call:  
**+1-608-829-8500**

Grass Disposable Deep Cup Electrodes continue a tradition of quality and innovation from the most trusted brand in EEG electrodes. Increased cup volume holds more paste for extended EEG recording.

- 50% more surface contact area and 20% greater cup depth than standard models
- Low impedance with both routine EEG and long term monitoring
- No cleaning and maintenance; use once and dispose
- Touchproof connector fit is optimized for easy headbox insertion and removal
- Packs of 10: ten unique lead wire colors
- Packs of 25: twelve unique lead wire colors (two each plus one extra white)




#### Grass Disposable Deep Cup Electrodes

Part Number	Cup Type	Lead Length	Quantity
023146	10mm Ag/AgCl	39" (1.0m)	10/pkg
023147	10mm Ag/AgCl	59" (1.5m)	10/pkg
023148	10mm Ag/AgCl	79" (2.0m)	10/pkg
023149	10mm Ag/AgCl	98" (2.5m)	10/pkg
027924	10mm Ag/AgCl	39" (1.0m)	25/pkg
027925	10mm Ag/AgCl	59" (1.5m)	25/pkg
027926	10mm Ag/AgCl	79" (2.0m)	25/pkg
027927	10mm Ag/AgCl	98" (2.5m)	25/pkg

**natus**®

# Grass® Disposable Standard Cup Electrodes



## Trusted, Reliable, Convenient

### New 25-Packs!

#### Materials

##### Electrode

Silver/Silver Chloride coated ABS plastic

##### Cup Dimensions

Cup diameter ..... 10mm

Cup height ..... 2.5mm

Interior cup diameter ..... 4.6mm

##### Lead Wire

Core ..... Tin plated Copper

Jacket ..... PVC

##### Touchproof connector, 1.5mm female (DIN 42 802)

Connector ..... Tin plated Brass

Insulation ..... PVC



##### RoHS compliant

To learn more about Natus products, contact your local distributor or sales representative.

US customers call:  
**1-800-356-0007**

International customers call:  
**+1-608-829-8500**

Grass Disposable Standard Cup Electrodes continue a tradition of quality and innovation from the most trusted brand in EEG electrodes. The standard cup electrodes are designed with a wide, flat rim for gentle contact with sensitive skin.

- Low impedance with both routine EEG and long term monitoring
- No cleaning and maintenance; use once and dispose
- Touchproof connector fit is optimized for easy headbox insertion and removal
- Packs of 10: ten unique lead wire colors 
- Packs of 25: twelve unique lead wire colors (two each plus one extra white) 

Grass Disposable Standard Cup Electrodes			
Part Number	Cup Type	Lead Length	Quantity
019-477200	10mm Ag/AgCl	39" (1.0m)	10/pkg
019-477300	10mm Ag/AgCl	59" (1.5m)	10/pkg
019-477400	10mm Ag/AgCl	79" (2.0m)	10/pkg
019-479200	10mm Ag/AgCl	98" (2.5m)	10/pkg
027920	10mm Ag/AgCl	39" (1.0m)	25/pkg
027921	10mm Ag/AgCl	59" (1.5m)	25/pkg
027922	10mm Ag/AgCl	79" (2.0m)	25/pkg
027923	10mm Ag/AgCl	98" (2.5m)	25/pkg



## Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. P/N 023765D

**natus**

Natus Medical Incorporated

natus.com