



Comprehensive, flexible and secure sleep diagnostic software

What's new in SleepWorks

- Windows® 10 Support
- New hardware platform
- Improved central administrator functions
- Automated observation notes
- EDF/EDF+ import and export

Natus SleepWorks is a comprehensive, flexible and cyber-secure sleep diagnostic solution that simplifies the setup, collection, review, analysis and reporting of comprehensive sleep studies. With advanced HIS/EMR integration and unmatched reliability, SleepWorks is a complete diagnostic package designed to meet the needs of private practice clinics, sleep clinics, hospitals and academic medical centers.

Enhanced workflow for sleep studies

- Ongoing study parameters at the forefront of view in a configurable Heads Up Display
- Configurable observation notes – automate tech notes with custom form fields and flexible reporting
- Analyze and review an ongoing studies, including video
- Capture PAP treatment and supplemental O₂ levels easily using customizable therapy types that simplify titration and reporting

Intuitive Clinical Experience

- Powerful event detection, trending and analysis options with fully synchronized HD video
- Record data directly to a server and access studies from any networked station
- Reduce size of archived files using advanced clipping and pruning features
- Use the Natus Database to quickly locate patient studies based on a variety of custom parameters



SleepWorks provides unmatched reliability and fault tolerant distributed networking solutions

Robust Security Features

- Meets HIPAA requirement via additional XLSecurity layer
- Define access to features, studies and programs based on logged-in user or role
- Automated power recovery
- Microsoft® BitLocker® drive encryption support

Site Management tools for enhanced productivity

- Enhanced security features with Active Directory integration
- Credentials cache system
- Automated site-wide software deployment with optional Cerebrum Enterprise Solution

A new cutting-edge hardware platform bringing together long standing reliability and advanced technology

Natus' new hardware platform is inspired by a rich heritage of expertise, drawing from our top sleep brands. The amplifiers consist of two main parts, a breakout box and a wall or cart mounted universal Natus Base Unit.

Breakout Box Features

- Small, rugged and lightweight amplifiers
- Powered sensor inputs
- Integrated pressure sensor
- Integrated pulse oximeter
- Up to 103 channels

Natus Base Unit Features

- Touchscreen display for real time interaction with ongoing studies
- Visual impedance checking and adding notes
- Perform bio-calibrations while in the patient room
- Universally compatible with several Natus EEG/PSG amplifier breakouts

Full range of powerful amplifiers



Standard sleep testing with **Embla SDx**



Advanced neuro-sleep monitoring with **Embla NDx**



Now integrated with **Embletta MPR Scalable Sleep System**

To learn more about Natus products, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. P/N 024435B

natus®

Natus Medical Incorporated

natus.com