







RIP* effort belt solutions sized to fit a wide range of patients

The Natus XactTrace family of respiratory effort belts provide a wide range of sizing options for in-home and in-lab sleep studies. XactTrace belts are reliable and provide consistently high quality data.

- 1 **Single Use Pre-sized Belts** provide ease of use, fast application and quality data
- 2 **Single Use Cut-To-Fit Belts** allow customized sizing to accommodate hard-to-fit patients
- 3 **Reusable Belts** provide quality data, are economical, and fit most patients

1 Single Use Pre-Sized Belts			
	XactTrace Single Use Pre-Sized Universal RIP Belts	Compatible with most PSG systems that accept a 1.5mm female touchproof connector	<ul style="list-style-type: none"> Powered by an external RIP processor Pediatric, Small, Medium and Large sizes Use only XactTrace Pre-Sized Locks/Snap Sensors
	XactTrace Single Use Pre-Sized Embla RIP Belts	Compatible with select Embla models that accept key connectors. See the Neurodiagnostic Supplies Catalog for compatibilities for each XactTrace belt.	<ul style="list-style-type: none"> Powered by Embla PSG System Pediatric, Small, Medium and Large sizes Use only XactTrace Pre-Sized Locks/Snap Sensors
2 Single Use Cut-To-Fit Belts			
	XactTrace Single Use Cut-To-Fit Universal RIP Belts	Compatible with most PSG systems that accept a 1.5mm female touchproof connector	<ul style="list-style-type: none"> Powered by an external RIP processor Customized sizing for each patient; ideal for pediatric and bariatric use Use only XactTrace Cut-To-Fit Locks/Snap Sensors
	XactTrace Single Use Cut-To-Fit Embla RIP Belts	Compatible with select Embla models that accept key connectors. See the Neurodiagnostic Supplies Catalog for compatibilities for each XactTrace belt.	<ul style="list-style-type: none"> Powered by Embla PSG System Customized sizing for each patient; ideal for pediatric and bariatric use Use only XactTrace Cut-To-Fit Locks/Snap Sensors Optional XactTrace Cable Extensions
3 Reusable Belts			
	XactTrace Reusable Universal RIP Belts	Compatible with most PSG systems that accept a 1.5mm female touchproof connector	<ul style="list-style-type: none"> Integrated battery power Adjustable Fits most patients
	XactTrace Reusable Embla RIP Belts	Compatible with select Embla models that accept key connectors. See the Neurodiagnostic Supplies Catalog for compatibilities for each XactTrace belt.	<ul style="list-style-type: none"> Powered by Embla PSG System Small, Medium and Large sizes

*Respiratory Inductance Plethysmography

Ordering Information

Additional XactTrace product information is available in the **Neurodiagnostic Supplies Catalog** which is available at natus.com

1 XactTrace Single Use Pre-Sized UNIVERSAL RIP Belts, Snap Sensors & RIP Processors

Part #	Type/Description
011516	Pre-Sized Single Use Belt, Pediatric
011517	Pre-Sized Single Use Belt, Small
011518	Pre-Sized Single Use Belt, Medium
011520	Pre-Sized Single Use Belt, Large
011521	Snap Sensor for Pre-Sized Single Use Belts, Abdomen, 26" (66cm) cable
011522	Snap Sensor for Pre-Sized Single Use Belts, Thorax, 16" (41cm) cable
1421048	RIP Processor Sensor, Abdomen, 78" (200cm) cable
1421047	RIP Processor Sensor, Thorax, 78" (200cm) cable

XactTrace Single Use Pre-Sized RIP Belts & Snap Sensors for EMBLA PSG Systems

Part #	Type/Description
011516	Pre-Sized Single Use Belt, Pediatric
011517	Pre-Sized Single Use Belt, Small
011518	Pre-Sized Single Use Belt, Medium
011520	Pre-Sized Single Use Belt, Large
011521	Snap Sensor for Pre-Sized Single Use Belts, Abdomen, 26" (66cm) cable
011522	Snap Sensor for Pre-Sized Single Use Belts, Thorax, 16" (41cm) cable
013224	Snap Sensor for Pre-Sized Single Use Belts, Abdomen, 78" (200cm) cable
013225	Snap Sensor for Pre-Sized Single Use Belts, Thorax, 78" (200cm) cable

2 XactTrace Single Use Cut-To-Fit UNIVERSAL RIP Belt Kits

Part #	Type/Description
1421052	Universal Single Use Cut-to-Fit Belt Kit - Complete*
1421050	Universal Single Use Cut-to-Fit Abdomen Kit*
1421051	Universal Single Use Cut-to-Fit Thorax Kit*
XactTrace Single Use Cut-To-Fit RIP Belts Components for EMBLA PSG Systems	
Part #	Type/Description
1421020	Lock for Cut-To-Fit Single Use Belt, ABD, 27.5" (70cm) cable
1421022	Lock for Cut-To-Fit Single Use Belt, Thorax, 15.7" (40cm) cable
1421021	Lock for Cut-To-Fit Single Use Belt, ABD, 27.5" (70cm) cable**
1421023	Lock for Cut-To-Fit Single Use Belt, Thorax, 15.7" (40cm) cable**
1421000	Cut-To-Fit Single Use belt material, 65.5' (20m)
1451070	Cable Extension, Abdomen, 78" (200cm), Optional
1451072	Cable Extension, Thorax, 78" (200cm), Optional

3 XactTrace Reusable UNIVERSAL RIP Belt Kits

Part #	Type/Description
1421091	Universal Kit - Complete
1421080	Universal Abdomen Kit (belt, sensor, cable)
1421081	Universal Thorax Kit (belt, sensor, cable)
XactTrace Reusable RIP Belts for EMBLA PSG Systems	
Part #	Type/Description
1421032	Small ABD, 22" (57cm) to 41" (104cm), 27" (70cm) cable
1421033	Small TRX, 22" (57cm) to 41" (104cm), 15.7" (40cm) cable
1421034	Medium ABD, 33" (83cm) to 56" (142cm), 27.5" (70cm) cable
1421035	Medium TRX, 33" (83cm) to 56" (142cm), 15.7" (40cm) cable
1421036	Large ABD, 47" (120cm) to 76" (194cm), 27.5" (70cm) cable
1421037	Large TRX, 47" (120cm) to 76" (194cm), 15.7" (40cm) cable
1421042	Small ABD, 22" (57cm) to 41" (104cm), 78" (200cm) cable
1421041	Small TRX, 22" (57cm) to 41" (104cm), 78" (200cm) cable
1421039	Medium ABD, 33" (83cm) to 56" (142cm), 78" (200cm) cable
1421038	Medium TRX, 33" (83cm) to 56" (142cm), 78" (200cm) cable
1421044	Large ABD, 47" (120cm) to 76" (194cm), 78" (200cm) cable
1421043	Large TRX, 47" (120cm) to 76" (194cm), 78" (200cm) cable

*Kits include locks, RIP processors and belt material

**For Embletta Gold only

Supplies

Convenient, complete, trusted

Natus supports the full spectrum of Sleep care, providing a complete portfolio of Sleep supplies for a seamless solution.

We also offer:

- Dedicated and knowledgeable customer service
- Streamlined order processing
- Convenient online ordering (US Customers only)

[Natus Medical Store – natusmedicalstore.com](http://NatusMedicalStore.com)

To learn more about Natus products, contact your local distributor or sales representative.

US customers call: **1-800-356-0007**

International customers call: **+1-608-829-8500**

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. **011236G**

natus

Natus Medical Incorporated

natus.com