



Compact EMG/NCS/EP Workstation



Dantec Keypoint Focus workstation



Dantec Keypoint Focus portable notebook system

The Dantec Keypoint Focus EMG system is based on nearly 50 years of experience in EMG system design. Quality, accuracy, and dependability – what you have come to expect from the Dantec brand.

Dantec Keypoint Focus is the next generation EMG system – delivering the quality and results you expect, at an economical price.

Its new compact design and intuitive, powerful functionality make the Dantec Keypoint Focus EMG system ideal for space conscious professionals in hospitals and private offices.

- All-in-One EMG/EP system with dedicated controls, integrated stimulators and high quality audio
- 3-, 4-, 6- or 8-channel amplifier options
- Advanced features for quantitative EMG, SF-EMG and tremor analysis
- On-line comparison to reference values

Handle difficult cases more quickly and efficiently

- Protocol-based testing
- Broad range of powerful analysis tools
- Automatic comparison with normative data
- The only system available with a complete MUP normative database
- Fast, automatic report generation

Design an EMG system that's perfect for you

- Select the configuration that works for you - portable notebook systems to fully featured workstations
- Purchase only the applications you need
- Completely control all measurements, data tables/calculations and displays
- Create customized reports


IT-friendly features

- Microsoft® SQL database
- Automated study data back-up/archiving functionality
- Complete Active Directory Services
- GDT/HL7 connectivity

Versatile EMG/NCS/EP Software




Multimodality NCS with reference values




Advanced H Reflex




Advanced Decrement test



Multi-MUP with reference values



Peak-triggering SF-EMG



EP with reference values

Service

At Natus, we strive for excellence in customer and technical service.

Here's how we can help:

- Accessible and effective Technical Support
- Definitive technical documentation and knowledgeable installation teams
- Replacement unit and spare part availability
- Extended warranty and service coverage programs
- Comprehensive, flexible customer training courses

Supplies

Convenient, complete, trusted

Natus supports the full spectrum of neurodiagnostic care, providing a complete portfolio of Sleep supplies for a seamless solution.

We also offer:

- Dedicated and knowledgeable customer service
- Streamlined order processing

To learn more about Natus products, contact your local distributor or sales representative.

International Customers Call: +1-608-829-8500

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. 007379F

natus

Natus Medical Incorporated

natus.com