



Quality. Comfort. Value.



Tapered, ultra-sharp,
3-bevel tip

0.38mm diameter

- Market pricing
- Excellent quality; ease of insertion; reliable signal
- Trusted source for neurodiagnostic instruments and IOM supplies
 - Precision cut, gently tapered, ultra-sharp lancet tip for effortless skin penetration assuring optimal patient comfort
 - High-tensile-strength, medical grade stainless steel needle for optimum strength and flexibility
 - Low profile needle-to-wire joint reduces pressure points on the patient's skin
 - Multi-colored lead wires with matching connectors for easy placement identification
 - Low noise and impedance level for clear, reliable recordings
 - 1.5mm insulated female connector, touchproof (DIN 42 802)
 - Individually packaged in sterile, pre-printed pouches for safety and easy identification

Technical Specifications

Part Number	Type	Needle Length/Diameter	Lead Wire Length	Quantity
019-476600	Stainless Steel, 6-color	13mm, 27g (0.38mm)	1m/39"	1/pouch, 24/box
019-476800	Stainless Steel, 6-color	13mm, 27g (0.38mm)	1.5m/59"	1/pouch, 24/box
019-477000	Stainless Steel, 6-color	13mm, 27g (0.38mm)	2.5m/96"	1/pouch, 24/box
019-476700	Stainless Steel, 12-color	13mm, 27g (0.38mm)	1m/39"	1/pouch, 24/box
019-476900	Stainless Steel, 12-color	13mm, 27g (0.38mm)	1.5m/59"	1/pouch, 24/box
019-477100	Stainless Steel, 12-color	13mm, 27g (0.38mm)	2.5m/96"	1/pouch, 24/box

4 each of 6 Lead Wire Colors: Red, Green, White, Black, Blue, Yellow



2 each of 12 Lead Wire Colors: Red, Green, White, Black, Blue, Yellow, Gray, Orange, Brown, Purple, Light Blue, Pink



Materials:

Needle Medical grade Stainless Steel
 Wire..... Copper
 Wire Insulation.....PVC
 Connector Type Touchproof, 1.5mm insulated female (DIN 42 802)
 Sterilization TypeEtO
 Shelf Life.....3 years for unopened pouches



Natus Ultra Subdermal touchproof connector

Supplies

Natus offers a full range of neurodiagnostic supplies, featuring:

- Dedicated and knowledgeable customer service
- Streamlined order processing
- Convenient online ordering (US Customers only)
Natus Medical Store – natusmedicalstore.com

To learn more about Natus products, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. **014742C**

natus®

Natus Medical Incorporated

natus.com