

Grass® Reusable "Tangle-free" Cup Electrodes



Durability, Ease of Use and Reliable Signals Time After Time



- Choice of gold, silver or silver/silver chloride 10mm cups
- Pediatric 6mm cup silver/silver chloride available
- No-lift contour adjustable cup neck for low profile (gold & silver cups only)
- Molded plastic 1.5mm female touchproof connector
- Low insertion force connector for easy headbox insertion and removal
- Six lead wire colors (two each): red, blue, green, white, brown, yellow

Materials

Electrode

Gold, silver or silver/silver chloride

Lead Wire

Core Kevlar and Silver plated Copper
 JacketTPU

Touchproof connector, 1.5mm female (DIN 42 802)

Connector Gold plated Brass
 Insulation..... PVC

Latex free / RoHS compliant

To learn more about Natus products, contact your local distributor or sales representative.


US customers call:
1-800-356-0007

International customers call:
+1-608-829-8500

Part Number	Cup Type	Lead Length	Quantity
019-477500	10mm Gold	39" (1.0m)	12/pkg
019-477600	10mm Gold	59" (1.5m)	12/pkg
019-478600	10mm Gold	79" (2.0m)	12/pkg
019-477700	10mm Gold	98" (2.5m)	12/pkg
019-477800	10mm Silver	39" (1.0m)	12/pkg
019-477900	10mm Silver	59" (1.5m)	12/pkg
019-478700	10mm Silver	79" (2.0m)	12/pkg
019-478800	10mm Silver	98" (2.5m)	12/pkg
019-478400	10mm Ag/AgCl	39" (1.0m)	12/pkg
019-478500	10mm Ag/AgCl	59" (1.5m)	12/pkg
019-478900	10mm Ag/AgCl	79" (2.0m)	12/pkg
019-479000	10mm Ag/AgCl	98" (2.5m)	12/pkg
019-479100	6mm Ag/AgCl Pediatric	59" (1.5m)	12/pkg



Grass® Disposable Cup Electrodes



Safety, Efficiency and Reliable Signals

Materials

Electrode

Silver/silver chloride

Lead Wire

Core Tin plated Copper

Jacket PVC

Touchproof connector, 1.5mm female (DIN 42 802)

Connector Tin plated Brass

Insulation PVC

Latex free / RoHS compliant

- Silver/silver chloride 10mm cups
- Single use eliminates the risk of cross-contamination
- No cleaning and maintenance; use once and dispose
- Flexible lead wires terminate in 1.5mm female touchproof connectors
- Five lead wire colors (two each): green, blue, red, yellow, white



Part Number	Cup Type	Lead Length	Quantity
019-477200	10mm Ag/AgCl	39" (1.0m)	10/pkg
019-477300	10mm Ag/AgCl	59" (1.5m)	10/pkg
019-477400	10mm Ag/AgCl	79" (2.0m)	10/pkg
019-479200	10mm Ag/AgCl	98" (2.5m)	10/pkg

Supplies

Natus offers a full range of neurodiagnostic supplies, featuring:

- Dedicated and knowledgeable customer service
- Streamlined order processing
- Convenient online ordering (US Customers only)
Natus Medical Store – natusmedicalstore.com

To learn more about Natus products, contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. **015468B**

natus®

Natus Medical Incorporated

natus.com