

# The Embletta® Multi Parameter Recording Family



## The scalable sleep system that grows with your needs

### ST Proxy

- 2 EEG
- 2 EOG
- Bipolar Chin EMG
- ECG
- 2 X Leg EMG



### ST+ Proxy

- 6 EEG
- 2 EOG
- Tri-polar Chin EMG
- ECG
- 2 X Leg EMG



### Embletta MPR PG

- Abdominal effort
- Thoracic effort
- Nasal Pressure Flow
- Snore
- SpO<sub>2</sub>
- Heart rate
- 1 x ExG
- Position
- Activity
- 1 x DC channel
- Microphone (optional)

### TX Proxy

- Differential Pressure
- 6 x DC (-5/+5V),
- Ambient Light
- Ability to record Flow/Pressure/Leak



- Scalable from Level IV to Full Level I PSG
- Multi-night recording for HST
- LCD screen for Patient Set-Up, Signal Quality Check and Viewing Traces
- Bio-Sensor Auto Start
- Software included
- Hi-Quality Microphone (optional)
- TX Proxy adds interfaces for DC devices

To learn more about Natus products, contact your local distributor or sales representative.

US Customers Call:  
1-800-356-0007

International Customers Call:  
+1-608-829-8500

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. 013165C

**natus**®

Natus Medical Incorporated

natus.com