



## A common software platform for EEG, LTM, ICU, Sleep and research studies

### What's new in NeuroWorks

- Advanced Trending
- Physiological Monitoring Integration
- Support for new LTM data management
- Improved display of annotations and navigation to regions of interest across a set of studies
- Windows® 10 compatibility
- Support for Quantum on-board memory and Nicolet Cortical Stimulator integrations
- Support for new cutting-edge Brain Monitor
- Frequency Analysis Tool
- Addition of Full HD (1080p) video

The Natus NeuroWorks platform simplifies the process of collecting, monitoring and managing data for routine EEG testing, long-term monitoring, ICU monitoring and Sleep studies. NeuroWorks systems are scalable to meet the needs of private practice clinics, hospitals and large teaching facilities.

#### Intuitive Clinical Experience

- Multiple Spike & Event analyzer options
- Fully synchronized HD video
- Monitoring options for home ambulatory studies
- Scalable remote access with Citrix provisioning
- Powerful qEEG trending

#### Time-Saving Administrator Functions

- Ease of setup using IP streaming cameras and amplifiers
- Powerful enterprise level SQL database
- Automated discovery of amplifiers and cameras using DHCP

#### Streamlined Workflow

- Easy data searching with a variety of customizable parameters
- Automated settings synchronization across networked NeuroWorks stations
- Convenient import and export with EDF/EDF+

## NeuroWorks configures to your network and security needs

### Robust Security Features

- HIPAA requirements via dedicated XLSecurity layer
- Role-based security access to features and studies
- Automated power recovery
- Microsoft® BitLocker™ drive encryption support

### Advanced IT Capabilities

- Bi-directional HL7 integration with HIS/EMR
- Citrix® and VMWare® Ready partner – the most powerful infrastructure to facilitate collaboration



## A cutting-edge, single solution for EEG, LTM, ICU and Sleep Studies

### Advanced software for clinical excellence

- EEG/Sleep combo software package
- User customizable personal workspaces
- Software Development Kit for direct data access for research
- Data-Stream Wizard with two data streams for efficient clinical and research workflow\*

### Reliable and flexible hardware options

- Patient wearable compact and lightweight amplifiers
- Range of EEG amplifiers with 32- to 256- channels, plus differential and DC channel options
- Wireless amplifier options for long-term monitoring
- Range of sampling rates from 200 Hz to 16 kHz
- Perform EEG recordings in DC mode\*\*
- Advanced Digital Trigger Input (up to 256 trigger values)\*\*

## A family of amplifiers to meet your needs



**NEW!** Natus Quantum



EMU40EX



**NEW!** Natus Brain Monitor



Nicolet® v32



Nicolet® v44



EEG32U



Trex HD

To learn more about Natus products contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**

\*Only available with Quantum

\*\*Only available with Quantum and Natus Brain Monitor

## Healthcare solutions with one thing in mind. You.

©2021 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. **013148 RevG**

# natus®

Natus Medical Incorporated

[natus.com](http://natus.com)