



REV	DESCRIPTION	CN	BY	APPR	DATE
B	UPDATE LAYOUT TO ALLOW LOCATION FOR UDI LABEL	21943	MPR	MS	22-JUN-16
C	1. CHANGED NATUS NEUROLOGY TO NATUS 2. UPDATED RX ONLY SYMBOL AND LEGAL MANUFACTURER DETAILS. 3. ADDED CAUTION, FOLLOW IFU, DO NOT USE IF PACKAGE IS DAMAGED, NOT MADE WITH NATURAL RUBBER LATEX, MANUFACTURING DATE, LOT NUMBER AND DISPOSAL SYMBOLS AND MEDICAL DEVICE TEXT. 4. ADDED TEMPERATURE LIMIT AND HUMIDITY LIMIT WITH ASSOCIATED SYMBOLS.	33941	TAHIR	SEE AGILE	SEE AGILE
D	UPDATED PRODUCT NAME, UDI INFORMATION AND LIFECYCLE PHASE CHANGED FROM PROTOTYPE TO PRODUCTION	35367	TAHIR	SEE AGILE	SEE AGILE

-- RoHS Compliant to 2011/65/EU* -- Latex Free -- BPA Free -- DEHP <0.1% by Weight

FRONT

BACK



LABEL NOTES:

1. TEXT AND GRAPHICS: INDELIBLE BLACK.
2. DO NOT PRINT RED BREAKING LINES AND RED TEXT.
3. DO NOT PRINT DIMENSION LINES.
4. *FINISHED PARTS TO COMPLY TO EU ROHS DIRECTIVE 2015/863
5. GLUE THE UDI BARCODE LABEL FROM 019615.
6. DATE OF MANUFACTURING WITH ASSOCIATED INFORMATION AND LOT NUMBER WITH ASSOCIATED INFORMATION ARE PRESENT ON THE LABEL 019615.

THIS DOCUMENT CONTAINS CONFIDENTIAL OR PROPRIETARY INFORMATION OF NATUS NEUROLOGY INCORPORATED. NEITHER THIS DOCUMENT NOR THE INFORMATION HEREIN IS TO BE REPRODUCED, DISTRIBUTED, USED OR DISCLOSED, IN PART OR IN WHOLE, EXCEPT AS SPECIFICALLY AUTHORIZED BY NATUS NEUROLOGY INCORPORATED.

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS. TOLERANCES ARE:		MATERIAL			
$.X \pm 0.1$ $.XX \pm 0.02$ $.XXX \pm 0.015$		MATERIAL #80			
$\angle \pm 1^\circ$ SURFACE FINISH 		FINISH		TITLE	
REMOVE ALL BURRS, BREAK SHARP EDGES DRAWING PER ANSI Y14.5 WELDING SYMBOLS PER AWS A2.0		GLOSS COATED		LABEL FOR GENUINE GRASS® SILVER DISC ELECTRODES - REUSABLE	
DO NOT SCALE DRAWING		DESIGN BY		SIZE	
		MATTHEW KELLER		DRAWING/PART NO.	
		DATE		REV	
		28-MAY-2014		B 014180 D	
		APPROVED BY		SCALE: NTS	
		SEE AGILE		E-FILE: 014180.ai	
		DATE		SHEET 1 OF 1	
		SEE AGILE			