



## Providing technical quality, clinical efficiency and patient comfort




- Coated stainless steel needle with ultra-sharp tip increases ease of insertion and enhances patient comfort
- Low-friction coating facilitates ease of repositioning
- Color-coded by needle size for easy identification
- Consistently low impedance levels for accurate recordings
- Reusable, flexible lead cable included with each box
- Pre-Sterilized



## Technical Specifications

Part Number	Length	Needle diameter	Recording area (mm <sup>2</sup> )	Hub color	Quantity
902-DMF25-S	1" (25mm)	0.36mm (28G)	0.28	Red	1/pouch - 48/box
902-DMG25-S	1" (25mm)	0.46mm (26G)	0.34	Yellow	1/pouch - 48/box
902-DMF37-S	1.5" (37mm)	0.36mm (28G)	0.28	Orange	1/pouch - 48/box
902-DMG37-S	1.5" (37mm)	0.46mm (26G)	0.34	Green	1/pouch - 48/box
902-DMG50-S	2" (50mm)	0.46mm (26G)	0.34	Blue	1/pouch - 48/box
902-DMG75-S	3" (75mm)	0.46mm (26G)	0.34	Violet	1/pouch - 24/box



## Reusable Lead Cable

Included with TECA Elite Disposable Monopolar Needle Electrodes

Part Number	Length	Description	Electrode End	Cable End	Quantity
902-RLC-TP	24" (61cm)	Unshielded Cable	DMN connector	1.5mm touchproof	6/pkg



### Intended use

TECA Elite Disposable Monopolar Needle Electrodes are intended for use with recording, monitoring and stimulation/recording equipment for the stimulation/recording of biopotential signals including electromyography (EMG) and nerve potential signals.

### Materials

#### Electrode

Needle ..... Stainless Steel  
 Needle Coating ..... Ultralon PTFE Green  
 Hub ..... Brass, Copper/Nickel plating  
 Hub Color Cover ..... Acrylonitrile-Butadiene-Styrene (ABS)  
 Protective Tube ..... Polyethylene

#### Packaging

Pouch ..... Tyvek/EVA  
 Box ..... Cardboard

### Sterilization

Type ..... Gamma Irradiation  
 Shelf Life ..... 3 years

### Environmental

LATEX FREE / RoHS Compliant  
 Storage Temperature ..... 0 to 132° F (-17.7 to 55° C)

## Supplies

Natus offers a full range of neurodiagnostic supplies, featuring:

- Dedicated and knowledgeable customer service
- Streamlined order processing
- Convenient online ordering (US Customers only)  
[Natus Medical Store – natusmedicalstore.com](http://Natus Medical Store – natusmedicalstore.com)

To learn more about Natus products,  
 contact your local distributor or sales representative.

US Customers Call: 1-800-356-0007

International Customers Call: +1-608-829-8500

## Healthcare solutions with one thing in mind. You.

©2021 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. 015827 RevC

# natus®

Natus Medical Incorporated

[natus.com](http://natus.com)