



Continuous brain monitoring wherever and whenever you need it – changing the standard of care in the ICU



Powerful neurological monitoring for clinical confidence

- Assist in diagnosis of patient's cerebral function
- Many Advanced Trending options to assist with quick identification of pathological signals for immediate intervention
- Physiological Monitoring Integration options
- Online observation can help the physician identify irreversible brain damage

Security and peace of mind for IT professionals and clinicians

- XLSecurity Layer with role-based security package and drive encryption
- Searchable and powerful enterprise level SQL server distributed database
- Bi-directional HL7 integration with HIS/EMR
- Automated software installation for large sites with Cerebrum Enterprise Solution

Natus NeuroWorks software

- Simple and intuitive interface designed specifically for intensive care workflow
- Cart-mounted for portability or mounted directly on a column or wall
- Leading-edge software with powerful trending options including qEEG
- True TCP/IP amplifiers and cameras with DHCP support plug and play anywhere on the network
- Data sharing in multi-site environments

Brain monitoring using Natus' renowned NeuroWorks system provides critical and immediate patient cerebral function to ICU staff and remote monitors for physician review.



Natus Brain Monitor*

- Up to 64 referential, up to 12 differential, 16 DC channels
- 4 kHz sampling rate
- On-board impedance check
- Touchscreen display on base unit with real time interaction with ongoing studies
- Visual impedance checking and adding notes
- Perform bio-calibrations from the patient bedside

*requires NeuroWorks 9.0 or higher



Natus EMU40EX

- 40 AC channels, 4 DC channels, SpO₂
- 1 kHz sampling rate
- Bluetooth wireless technology with on-board memory for live monitoring
- Lightweight, patient wearable breakout box



Natus cybersecurity program

Proactively addressing cybersecurity as an integral part of our design process

- Secure handling of confidential information generated or maintained by our hardware and software
- Continuous evaluation of our products against identified threats and vulnerabilities in medical technology
- Monitoring of feedback channels to manage security events in the field

Natus support

At Natus, we strive for excellence in customer and technical service.

Here's how we can help:

- Accessible and effective Technical Support
- Definitive technical documentation and knowledgeable installation teams
- Optional extended warranty and service coverage programs
- Comprehensive, flexible customer training courses

EEG and ICU supplies solutions

Convenient, complete, trusted

Natus supports the full spectrum of EEG and ICU care, providing a complete portfolio of supplies for a seamless solution.

- Dedicated and knowledgeable customer support
- Streamlined order processing
- Convenient online ordering (US Customers only)
[Natus Medical Store – natusmedicalstore.com](https://natusmedicalstore.com)

To learn more about Natus products contact your local distributor or sales representative.

US Customers Call: **1-800-356-0007**

International Customers Call: **+1-608-829-8500**

Healthcare solutions with one thing in mind. You.

©2021 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. 024992 Rev03

natus®

Natus Medical Incorporated

natus.com