

REVISION HISTORY

CHG	EO	DESCRIPTION	DATE
A	ECO13776	Converted from D40274. Updated artwork revision and packaging requirements. Changed any Integra references to Natus.	25-JUN-2018
AA	ECO#31642	Change CE 0086 to CE 2729	12-Nov-2019
AB	ECO#33098	Change CE 0086 to CE 2797	11-Feb-2020
AC	ECO#36597	Remove CE Mark and EC REP	19-Mar-2021

1. DESCRIPTION:

Multi-Language DFU (English, French, German, Italian, Spanish, and Dutch), 110-4BT

2. MANUFACTURER:

A RR Donnelley

3. MANUFACTURER PART NUMBER:

A 40274

4. REFERENCE DOCUMENTS:

Q00035 Incoming Quality Assurance Operating Procedure
 Q00020 Receiving Inspection Procedure for Labeling Components From Outside Vendors
 A40274SD Rev. GD

5. SUPPLIER/DISTRIBUTOR:

B See approved supplier list

6. SUPPLIER/DISTRIBUTOR CATALOG NUMBER:

B See approved supplier list

7. MATERIAL SPECIFICATION:

80 lb. Matte Text or equivalent
 40-page booklet
 Double sided text including cover and back page

8. SHELF LIFE REQUIREMENTS:

N/A

9. VISUAL REQUIREMENTS:

G All printing to be 100% black tone on white background
 Inspect per Q00020

10. CLEANLINESS REQUIREMENTS:

E DFU shall be free of smears, dots and lines.

11. PACKAGING/LABELING REQUIREMENTS:

SOURCE CONTROL DRAWING



TITLE **DFU, Multi-Language, 110-4BT**

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF NATUS AND IS SUBJECT TO THE CONDITIONS THAT THE INFORMATION (A) BE RETAINED IN CONFIDENCE AND (B) NOT BE REPRODUCED OR COPIED IN WHOLE OR IN PART EXCEPT UNDER AN EXPRESS WRITTEN AGREEMENT WITH NATUS.

DWG NO **40274SD**

REV **AC**

ERP STOCK NO **40274SD**

FILE NAME N/A

SHEET **1 of 2**

D Single shrink wrapped (tight) Each box shall be labeled, at a minimum with Natus name, P/N, PO and Quantity.

12. STORAGE REQUIREMENTS:

May store in a warehouse environment. Store at ambient temperature and humidity, keep dry and clean.

13. SUPPLIER DATA COLLECTION FORM:

C Quality Assurance Report required with each shipment.

14. DIMENSIONAL:
(Finished Size)

H Length: 8.50"± 0.10"

Width: 5.50"± 0.10"

15. INSTRUCTIONS:

Sample booklet to be given to printer with artwork and drawing.
Finished booklet to be folded in center and stapled in 2 places at the spine.

SOURCE CONTROL DRAWING

DWG NO	40274SD	REV	AC
TITLE	DFU, Multi-Language, 110-4BT	SHEET	2 of 2