

This product is being supplied with an electronic copy of the Instructions for Use.

A copy the Instructions for Use in PDF format can be accessed by either locating it on the accompanying disk/USB Key or by contacting Natus Technical Support.

The files can be printed, saved, or searched using Adobe Reader. A copy of Adobe Reader has been included on the accompanying disk/USB Key or can be downloaded directly from Adobe Systems (www.adobe.com).

Description of Symbols that may appear on your system.

Symbol	Description	Symbol	Description
	ATTENTION: Consult Accompanying Documents		Dangerous Voltage
	Consult Operating Instructions. Failure to follow operating instructions could place the patient or operator at risk. Image on blue background.		Alternating Current
	Consult Accompanying Documents		Direct Current
	European Authorized Representative		Power On
	Manufacturer		Power Off
	Pushing Prohibited		EU only: Do Not Dispose as Unsorted Municipal Waste
	Protective Earth (Ground)		CE Mark
	Type BF Equipment (Isolated Inputs)		Class II Equipment (non-grounded enclosure)
	ESD Sensitive or Static Sensitive		DC Inputs 13 to 16 Breakout Cable Connector
	Base to Breakout Cable		Patient Event Button Connector
	Breakout to Breakout Cable Connector		Nonin Pulse Oximeter Sensor SpO2 Connector / Not for Continuous Monitoring
	Cortical Stimulator Adapter Connector		Equipotentiality
	Breakout Box Status LED		CAUTION: USA Federal law restricts this device to sale or on the order of a licensed Neurology practitioner.
	Impact hazard		Tipping hazard
	Minimum/maximum monitor weights		Cart transport
	Pinch hazard		Non-waterproof Device
	RF Equipment for Non-ionizing Radiation		Stored energy hazard
	Wireless Bluetooth Identification		
	ETL Certifications		Medical Device

To request a hard copy of the Instruction for Use, contact Technical Support via phone or email, or contact your local representative. Natus Technical Support is available 24 hours a day, 7 days a week.



Natus Medical Incorporated
 DBA Excel-Tech Ltd. (XLTEK)
 2568 Bristol Circle Oakville,
 Ontario, L6H 5S1 Canada
 Tel: 1-613-254-8877
 Toll Free (US & Canada): 1-888-662-7632
 Technical Support Email: Ottawa.TechnicalSupport@natus.com
 Website: www.natus.com



EUROPEAN AUTHORIZED REPRESENTATIVE
Natus Manufacturing Limited
 IDA Business Park, Gort,
 Co. Galway, Ireland

