



Streaming newborn videos for families

NICVIEW® 2 Live Video Streaming

Technical Specifications

Camera

Camera Body

Dimensions (HxLxD)	68 x 130 x 33 mm (2.7 x 5.1 x 1.3 in)
Weight	10 oz

Display

Type	Color, TFT, touchscreen
Dimensions	3.5"
Resolution	320x480 resolution

Operating Environment

Temperature Range	0-35° C (32-95° F)
Humidity Range	20-95% rel., non-condensing

Power Supply

Micro USB Power Supply	
Input Voltage/Range	100-240V AC, 50-60 Hz
Output Voltage	5.4V DC

Camera Lens

Effective Pixels	8 megapixels
Minimum Illumination	1 lx
Video Output	Up to 30 FPS at 1280x720p

Network

Ethernet	10/100 BaseT
Wireless Frequency	2.4 GHz / 5 GHz
Wireless Data Transfer	802.11.b/g/n/ac wireless
Wireless Security	WPA/WPA2-PSK (TKIP-AES)
Required Protocol	IPv4 - HTTPS
Suggested Protocol	IPv4 - UDP 3478 for highest quality video
Average bandwidth consumption during viewing	500 KBps
Average bandwidth consumption when not viewing	1.5 KBps

Client Access – Supported Browsers

On Microsoft Desktop OS

- Google Chrome (latest release version)
- Firefox (latest release version)
- Internet Explorer 10 and 11 (with the OpenTok Plugin for IE)
- Opera (latest release version)
- Safari (latest release version)

On macOS and iOS

- Safari (latest release version)
- iOS NICVIEW app available on the Apple App Store

On Android

- Google Chrome for Android
- Firefox for Android
- Android NICVIEW app available on the Google Play Store


IMPORTANT – SAFETY PRECAUTION:

Carefully read the instructions for use before using this product.

For Technical Support, please contact us at <https://www.nicview.net/help.aspx>

Mounting Arm

Mounting Arm	31" long
Arm Weight	3.2 lb



Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. P/N 026839B

natus

Natus Medical Incorporated

[natus.com](https://www.natus.com)