



Providing exceptional eye protection for neonatal phototherapy



"We liked the Biliband (eye protector) better than the goggles our first preemie had to use! The goggles didn't stay on nearly as well as the Biliband."

– Jason W., Greenwood, IN

"We like the new Bilibands because they stay in position better than our previous products and they are gentle on the baby's skin."

– Children's Hospital, New Orleans, LA

Blocks harmful light to baby's eyes

- Provides protection from phototherapy light
- Special eye-pad material blocks light to the eyes, while headband design allows phototherapy light to reach baby's head

Provides a perfect fit

- Remains securely in place with unique Y-shaped design
- Conforms to any head shape with easy-to-use, adjustable tabs

Safe and comfortable

- Biliband material minimizes the risk of skin irritation and is not made with natural rubber latex
- Feels soft and comfortable around baby's head and face

Available in 3 sizes:

Ordering Information

Size	Head Circumference	Part #	Packaging
Regular	30 - 38 cm	900642	20 units/pk
Premature	24 - 33 cm	900643	20 units/pk
Micro	20 - 28 cm	900644	20 units/pk

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. 051262J

natus

Natus Medical Incorporated

natus.com