



## Comprehensive care from screening to early intervention

Neometrics iCMS provides a robust data management system that is integral to a state's ability to successfully coordinate the time-sensitive follow-up activities needed to ensure a child gets from screening to diagnosis and into early intervention.

### Key features

- Neometrics iCMS software integrates case management for multiple tracking and monitoring programs including birth defects, hearing screening, lead screening, critical congenital heart disease (CCHD) and Dried Blood Spot screening
- iCMS software tracks infants from initial entry into the system through receipt of services to case closure. The system can take in referral results, audiologist information, caregiver notes, appointment information and various result data
- iCMS software is able to automatically assign actions based on your algorithms and then can configure in the software and store. The software can base assignment of events on result values, dates, risk factors and can associate demographic influences into the final event to be assigned. The system has the ability to automatically notify caregivers by email of an event based on an incoming result or a manually entered event
- The system can provide long-term follow-up of patients after they have been put into treatment and track performance measures
- iCMS software currently utilizes HL7 for bidirectional communications with hospital information and medical records system. HL7 is a standard for data transmission that allows for demographic and result data to be transmitted utilizing the same standard

## Service highlights

- Natus is the leader in newborn hearing screening and newborn screening data management in the U.S.
- Neometrics has provided newborn screening solutions to approximately twenty (20) states over the past thirty (30) years
- Many of Neometrics' clients have been retained in excess of 15 years with some over 25 years

## Support & partnership

- Neonatal Care Academy – the online resource for education – [neonatalcareacademy.com](http://neonatalcareacademy.com)
  - Clinical white paper – Data management software and services: The impact of successful tracking and follow-up
  - Clinical updates
  - Clinical tool downloads



Join the Neonatal Care Academy for clinical resources, learning opportunities and more!  
[neonatalcareacademy.com](http://neonatalcareacademy.com)

## Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. 025892-001B

# natus®

Natus Medical Incorporated  
[natus.com](http://natus.com)