

# OraSwab™


## foam tip applicator

### Intended use:

General cleaning and removal of bacteria in the oral cavity.


### Instructions for use:


Using a moist swab, gently clean oral cavity to remove bacteria. Use hospital protocols for oral care.


 **Warning!** Avoid vigorous cleaning. Foam swab tip may separate from handle. Do not allow infant to suck on swab.


### SAFETY SYMBOLS USED:


#### Symbol Meaning

 Consult accompany documentation

 Do not reuse, single patient use only


 Batch code / lot number

 Expiration date / Use by date

 Caution and Warnings

**OraSwab 6mm and 9mm wide Foam Tip Applicators Tip is 20mm long with 155mm handle. Not made with natural rubber latex.**

Distributed by:  
Natus Medical Incorporated  
5900 First Avenue South  
Seattle, WA 98108 USA  
+1-650-802-0400  
1-800-303-0306 (US/CAN)  
www.natus.com

 Manufactured by:  
Puritan Medical Products Company LLC  
31 School Street, Guilford, ME 04443 USA

 EMERGO EUROPE  
Prinsessegracht 20, 2514 AP, The Hague  
The Netherlands



**natus**  
newborn care

021758

