



The ultimate standard for accurate weight



Smart Scale model 65



Smart Scale model 60

Olympic Smart Scale Models 60 and 65 combine remarkable accuracy with simple operation. To weigh a baby, simply place the child in the tray and the scale will automatically compute weight to the nearest gram – even if the baby is moving.

With the Smart Scale, changing numbers do not have to be tracked – the final weight is automatically locked on the display. One can continue caring for the baby and chart the weight when convenient. Even if the baby is removed from the scale, weight remains locked on the display.

- Weigh active babies quickly and accurately
- Four-sided easy lift trays are safer for babies while weighing
- Available as table-top or roll-around models

Accurate measurements

- Automatically computes weight to the nearest gram
- Automatically zeros itself before each weighing for accurate reproducible weight measurements
- Zero (tare) button allows extra weight, such as blankets or toys, to be cancelled out to accurately weigh only the baby

Easy to use

- **Fully automatic** after placing baby on the scale, weight is automatically computed and locked on the display
- **Automatic weight averaging** when the scale detects motion, it automatically averages and computes within seconds and locks the final weight on the display
- **Memory recall button** allows display of the last averaged weight, even if the scale has been turned off
- **Four-sided easy lift tray (Smart Scale only)** – safer for baby and hand ports make it easier for lifting baby in and out
- Smooth contours, molded shape and edgeless tray enable easy cleaning

Efficient

- Weigh babies within seconds, even when they are moving
- Weighs on rechargeable batteries or connected to an electrical outlet while batteries are charging
- Automatic shut down feature conserves battery life
- Available as either a table-top style or as a roll-around model on a sturdy cart

Roll-Around Smart Scale®

Model 65

- Patient is protected with scale securely anchored to the sturdy roll-around base
- Padded crash-bar handles at both ends protects the scale when moved around
- Lower height makes lifting infants in and out of the scale safer for the baby, while reducing potential back strain for the clinician
- Four fully-swiveling, ball-bearing casters enable easy maneuvering
- Supplies can be stored in large enclosed cabinet with transparent doors



Contact your local Natus representative
call 1-800-356-0007

*Product not available in all markets. Contact your local representative for additional information.

Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. 011500 RevB

natus®

Natus Medical Incorporated

natus.com