

Compact warmers for nurseries,
birthing rooms, emergency
departments, operating rooms
and recovery rooms

Safe and convenient warming of
blankets and external solutions

Maintain a supply of warm blankets
near your patients

Compact size and silent operation
are ideal for nursery environment



Warmette Model 10



Warmette Model 20

Space efficient

- Compact exterior with large interior capacity
- Sized to position on counter, mount to wall, or on sturdy floorstand
- Use side-by-side for additional storage or for dual temperatures

Safe and effective

- Lockable temperature setting
- Optional feature allows temperature setting knob to be locked at five temperature set points (100°F, 110°F, 120°F, 130°F, and 140°F)
- Radiant heating element provides safe, uniform, controllable temperature
- Silent operation is ideal for nursery environment
- Rounded shape has no sharp corners

User-friendly controls

- Simply dial the desired temperature

Signal light

- Indicates when chamber is above set temperature

Transparent door

- View blanket or external solution inventory at a glance

No special wiring needed

- Simply plug in and start using

Convenient

- Designed to position near patients to save time and steps

Flexible

- May easily be moved to accommodate future departmental changes

Features and Specifications*

Dimensions

Model 10: 16.5"H x 20"W x 16"D (41.9 x 50.8 x 40.6 cm)

Model 20: 23"H x 26"W x 24"D (57.8 x 66.5 x 61.3 cm)

Internal capacity

Model 10: 1.2 cubic ft. (0.03 m³)

Model 20: 4.0 cubic ft. (0.12 m³)

Construction

Double walled, painted aluminum exterior

Heating

Radiant heating element encapsulated in silicon – silent operation

Cleaning

Wipe down exterior and interior surfaces with ordinary liquid disinfectants

Adjustable temperature range

75°-150°F (24°-66°C)

Controls and signals

- Power switch
- Set temperature control
- Heater on
- Over set temperature lights – Indicate the internal temperature of the heating chamber is higher than set temperature

Electrical safety

All models CSA C/US certified to CSA C22.2 1010-1 and UL 61010A-1

Electrical Requirements

No special wiring required

Model 10: 120V~, 60Hz, 1.5A; or 220V~, 50Hz, 0.8A; 6' Power Cord

Model 20: 120V~, 60Hz, 3.5A; or 220V~, 50Hz, 1.9A; 6' Power Cord

Ordering Information

Olympic Warmette Model 10

For use on countertop, shelf, or table

CAT. NO. 56910

Wall Bracket for Warmette Model 10 (optional)

CAT. NO. 56917

Floor Stand for Warmette Model 10 (optional)

CAT. NO. 56916

Olympic Warmette Model 20

For use on countertop, shelf, or table

CAT. NO. 56920

Wall Bracket for Warmette Model 20 (optional)

CAT. NO. 56927

Floor Stand for Warmette Model 20 (optional)

CAT. NO. 56926

***Note:** Warmettes should be used to warm blankets and external-use solutions only

natus

Natus Medical Incorporated

1501 Industrial Road
San Carlos, CA 94070 USA
1-800-303-0306
+1-650-802-0400
www.natus.com

© 2013 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.

For Information or Ordering

Call 1-800-303-0306

(toll-free in US/Canada)