



Broaden Your Knowledge & Success at the Natus NERVE Center®

At Natus, we know education shapes the practice of many medical device-driven programs, from screening to diagnosis. Opportunities to improve health care have expanded, and so have our educational offerings.

Whether it's a policy/procedure, self-paced learning tutorial, product quick guide, or service manual, visit the Natus Electronic Resource Venue for Education, the NERVE Center, to satisfy your 'impulse' for learning! Our mission is to provide clinicians and biomedical engineers with advanced support and training materials needed for successful product use.

eSeminars

With the goal of delivering relevant education in a method that is easy and convenient, Natus offers monthly LIVE eSeminars focusing on cutting edge science and technology topics.

These interactive programs are provided by guest luminary medical speakers and in-house clinical subject matter experts, and are designed to contribute to your professional growth.

Can't attend the LIVE eSeminar? No need to worry. We've recorded the presentation for your 'ON DEMAND' viewing convenience. Watch and listen to a variety of archived eSeminars, at your pace, on your schedule!

Check out the Natus NERVE Center to register for LIVE online events or to view archived events.

Resource Library

Whether you're looking for a sample policy and procedure, parent education, or a Quick Start Reference Guide, the Natus NERVE Center Resource Library provides numerous reference pieces to meet your educational or procedural needs.

Take advantage of these tools to train new staff, expand staff development, and/or complete annual staff competencies. All items are cross-referenced by clinical area and training format to make your search a "synapse"!

Training Programs

Are you new to your field of practice? Or, just new to Natus technology? Browse our menu of learning paths where you'll find all of the relevant learning tools particular to your area of care. Participate in as many programs as you need, all at your own pace.

Natus is dedicated to providing education that enhances the technology solutions provided by our products. Broaden your knowledge by visiting the Natus NERVE Center at www.natus.com/nervecenter and find the educational tools you need to succeed!



Electrodes02

Skin Prep, Adhesives & Application.....04


Ear Tips & Cables06

Miscellaneous Supplies09

Natus, the leader in the hearing assessment market, offers integrated solutions for hearing screening and diagnostics. Nurses, audiologists, otolaryngologists and trained technicians use our systems for objective screening and diagnosis of hearing loss.

We offer a full range of hearing supplies and accessories for our Bio-logic hearing screening and diagnostic equipment.





Part #: 101001 - 101007



Part #: 101203 - 101204



Part #: 301003 - 301004

REUSABLE ELECTRODES

Natus Reusable EEG/EP 10mm Disc Electrodes

"High-dome" electrodes are available in different electrode materials: silver, silver/silver-chloride and gold as well as as two cable types: silicone (tangle-free) and Teflon (light-weight). These multi-colored leads end in a 1.5 mm molded connector.

Disc Electrodes with Silicone "Tangle-Free" Lead Wires

Part #	Item	Lead Length	Qty
101001	Silver Disc Electrode	1 m	10/pkg
101004	Silver-Chloride Disc Electrode	1 m	10/pkg
101007	Gold Disc electrode	1 m	10/pkg

Disc Electrodes with Teflon "Light-Weight" Lead Wires

Part #	Item	Lead Length	Qty
101001-T	Silver Disc Electrode	1 m	10/pkg
101007-T	Gold Disc Electrode	1 m	10/pkg

Reusable Ear Clip Electrodes

Two disc type electrodes mounted on a spring clip with an attached 1 meter lead (1.5 mm connector). Available in three metal types: silver, silver chloride and gold; and two disc sizes 6mm (pediatric) and 10mm (adult).

Part #	Item	Qty
101203	Silver Ear Clip, Ped, 6 mm	1/pair
101204	Silver Ear Clip, Adult, 10 mm	1/pair
101213	Silver-Chloride Ear Clip, Ped, 6 mm	1/pair
101214	Silver-Chloride Ear Clip, Adult, 10 mm	1/pair
101223	Gold Ear Clip, Pediatric, Ped, 6 mm	1/pair
101224	Gold Ear Clip, Adult, 10 mm	1/pair

Jumper Cables

The electrode jumper cables allow the same electrode to be plugged into two different inputs on the reamp. Individually packaged (1.5 mm connector).

Part #	Item	Qty
301003	Reusable Jumper Electrode - Red	1
301004	Reusable Jumper Electrode - Black	1

REUSABLE ELECTRODE LEADS

Lead Wires for use with Snap Disposable Electrodes

Part #	Item	Lead Length	Qty
301601	Snap Lead - (R/WT/BL/G/BK)	1 m	5/pkg
301663	Pinch Clip Leads - (R/WT/BL) (for snap electrodes)	18 in	3/pkg
301666	Pinch Clip Leads - (R/WT/BL/G/BK/Y) (for snap electrodes)	18 in	6/pkg

For use with Tab Style Disposable Electrodes

Part #	Item	Lead Length	Qty
301633	Alligator Clip Leads - (R/WT/BL)	18 in	3/pkg
301636	Alligator Clip Leads - (R/G/BK)	1 m	3/pkg
301635	Alligator Clip Leads - (R/WT/BL/G/BK)	18 in	5/pkg

DISPOSABLE ELECTRODES

Disposable Adhesive Electrodes

Custom-designed electrodes provide optimum performance and comfort to meet the needs of a varied patient population. Three-pack foil pouches keep them fresh.

SINGLE PATIENT USE ONLY.

Part #	Item	Qty
101600	Snap electrodes	60/pkg
101603	Snap Electrodes	300/bx
101604	Square Snap Electrodes - aggressive adhesion (not for use on infants)	30/pkg
101605	Square Snap electrodes - aggressive adhesion (not for use on infants)	300/pkg
101661	Tab Electrodes	60/pkg

Subdermal Needle Electrodes - Disposable

Rhythmlink ultra-sharp 13mm 27-gage high tensile strength, vacuum melt, medical grade stainless steel needle electrodes, 1m lead wire.

6 color box contains 4 of each: white, red, green, dark blue, black and yellow lead wires

Part #	Item	Qty
102516	Subdermal Needle Electrodes - 6 color	24/bx

Speciality Electrodes - Tiptrodes Gold-Foil Foam Tips

Disposable electrode combined with an insert earphone foam tip used to record evoked potentials from the ear canal.

Part #	Item	Qty
101111	Tiptrode Gold-Foil Foam Tips - Pediatric	20/pkg
101113	Tiptrode Gold-Foil Foam Tip - Adult	20/pkg
301116	Tiptrode Electrode Cable/Acoustic Tube - (R/BL)	1/pair



Part #: 301601



Part #: 301633 - 301666



Part #: 301633 - 301635



Part #: 101604 - 101605



Part #: 101661



Part #: 102516



Part #: 101111, 101113, 301116

ELECTRODE SKIN PREP, GEL, PASTE & ADHESIVE PRODUCTS



Part #: 102566N

NuPrep® Skin Preparation

Topical skin abrasive. For low impedance. Mild abrasive for patient comfort. Chloride and acetone free.

Part #	Item	Qty
102566N	NuPrep Skin Prep - 4 oz tube	3/box
102566N-6	NuPrep skin Prep - 25 gram tube	6/box



Part #: 109083

Patient Prep Tape

Dry abrasive skin preparation. For low impedance.

Part #	Item	Qty
109083	Patient Prep Tape - Dispenser Roll - 196" roll	1



Part #: 109020

Alcohol Pads

Alcohol prep pads for antiseptic skin preparation prior to electrode application.

Part #	Item	Qty
109020	Alcohol Prep Pads - 200/box	1



Part #: 102463N

Ten20® Conductive Paste

Non-drying, water soluble conductive paste. Recommended for EP and EEG disc electrodes and also used under collodion electrodes.

Part #	Item	Qty
102463N	Ten-20 conductive Paste - 4 oz. tube	3/box



Part #: 109323 Part #: 109328

Lectron® II Gel

Chloride free, hypoallergenic and water soluble.

Part #	Item	Qty
109323	Lectron II Conductivity Gel - 3 oz. tube	1
109328	Lectron II Conductivity Gel - 8.45 oz. tube	1

3M™ SURGICAL TAPE

Micropore™ Surgical Tape

Micropore surgical tape is ideal for sensitive skin and elderly patients. It is an excellent surgical tape for repeated applications because it's highly permeable material allows moisture to evaporate permitting the skin to breathe. Miropore is easy to remove and it leaves minimal adhesive residue. Not made with natural rubber latex. Hypoallergenic.

Part #	Item	Qty
109401	3M Micropore Surgical Tape 1" (25 mm)	12 rolls/box



Part #: 109401

Transpore™ Surgical Tape

Transpore surgical tape is a transparent plastic tape perforated with a polyethylene film. This tape can be torn bi-directionally into very small strips even while wearing gloves. Strong adhesion keeps this tape in place. Not made with natural rubber latex. Hypoallergenic.

Part #	Item	Qty
109403	3M Transpore Surgical Tape 1" (25 mm)	12 rolls/box



Part #: 109403

Durapore™ Surgical Tape

Durapore is a "silk like" surgical tape that is an ideal general purpose tape for securing regular and bulky dressings, tubing and immobilization because of its high strength, high adhesion and durability.

Part #	Item	Qty
109405	3M Durapore Surgical Tape 1" (25 mm)	12 rolls/box



Part #: 109405

APPLICATION SUPPLIES

Tape Measure

The vinyl tape measure is marked in both centimeters and inches with automatic retraction at the push of a button. Measures up to 60" or 150 cm.

Part #	Item	Qty
109001	Tape Measure	1



Part #: 109001

Wax Marker

Great for skin marking. Water soluble - easy cleanup. Sold individually.

Part #	Item	Qty
109012	Wax Marking Pencil - Red	1



Part #: 109012

Ink Marker

The skin marker uses a brilliant blue ink visible on all skin tones.

Part #	Item	Qty
109013	Skin Marker	1



Part #: 109013

APPLICATION SUPPLIES



Part #: 109030

Cotton Tip Applicator

Cotton tipped applicator sticks (6" length). Non-sterile.

Part #	Item	Qty
109030	Cotton tip applicator	1000/bx



Part #: 109040

Gauze Sponge

Ideal for prepping, scrubbing and general purpose cleaning. Non-sterile, 4 ply.

Part #	Item	Qty
109040	Gauze Sponge, 2" x 2"	200/bx



Part #: 109051

Sani-Cloth®

These disinfectant towelettes are premeasured and kill microorganisms to reduce the risk of cross contamination. Gentle on hands. 135 wipes, 6" x 7.5" per container.

Part #	Item	Qty
109051	Sani-cloth Plus disposable Wipes (135 wipes per container)	1

OAE PROBE TIPS & ACCESSORIES



Part #: 203200 - 203219

Bio-logic® disposable OAE probe tips can be used for OAE testing with Bio-logic AuDX®, AuDX® Pro, Scout® and ABAer® systems. They can also be used for ABR testing with the 580-SINSER-012 insert earphones using the infant tip adapter (#202240). OAE probe tips come in a variety of sizes from preemie to adult to fit most ears and are conveniently packaged in 20/bag or 100/bag. Not made with natural rubber latex. *Single patient use only.*

Part #	Item	Outer Diameter	Length	Qty
203200	OAE Jumbo Tip	17.91 mm	12.70 mm	20/bag
203210	OAE Jumbo Tip	17.91 mm	12.70 mm	100/bag
203202	OAE Adult Foam Tip	13.85 mm	12.70 mm	20/bag
203212	OAE Adult Foam Tip	13.85 mm	12.70 mm	100/bag
203201	OAE Pediatric Foam Tip	9.50 mm	12.70 mm	20/bag
203211	OAE Pediatric Foam Tip	9.50 mm	12.70 mm	100/bag
203333	OAE 6 mm Blue Tip	6.0 mm	16.51 mm	20/bag
203353	OAE 6 mm Blue Tip	6.0 mm	16.51 mm	100/bag
203208	OAE Tree Tip	3.60 mm	15.10 mm	20/bag
203218	OAE Tree Tip	3.60 mm	15.10 mm	100/bag
203209	OAE Preemie Tree Tip	3.60 mm	15.10 mm	20/bag
203219	OAE Preemie Tree Tip	3.60 mm	15.10 mm	100/bag



Part #: 203000 & 580-PRBNZL

OAE Probe and Accessories

Part #	Item	Qty
203000	OAE Probe Test Cavity	5/pkg
580-PRBNZL	OAE Probe Body Nozzle	5/pkg
580-PROAE3	OAE Ear Probe Assembly	1



Part #: 580-PROAE3

INSERT EARPHONE DISPOSABLE TIPS & SUPPLIES

Part #	Item	Qty
202121	Insert Earphone Foam Tips - Pediatric	50/bag
202122	Insert Earphone Foam Tips - Adult	50/bag
202123	Insert Earphone Foam Tips - Jumbo	24/bag

NOTE: Compatible with 580-SINSER-012 Insert Earphone. Used with Bio-logic Navigator Pro® AEP and MASTER II systems.

Part #	Item	Qty
202240	Insert Earphone Infant Tips Adapter w/Acoustic Tubes *202240 shown with OAE probe tips	1/pair
202250	Insert Earphone Tubing Replacement	1/pair
202260	Insert Earphone Adapter tips with out Tubes	10/pkg

NOTE: Compatible with 580-SINSER-012 Insert Earphone. Used with Bio-logic Navigator Pro AEP and MASTER II systems.

Part #	Item	Qty
206920	Halo Ear Muffin Adapter Tubes	1/pair
206202	Halo Ear Muffin with Tab Electrodes	20 sets/box

NOTE: Compatible with 580-SINABR-008 and 580-SINSER-012 Insert Earphones. Used with Bio-logic ABAer and MASTER II systems.



Part #: 202121 - 202123



Part #: 202240*



Part #: 202250

Part #: 206920

Part #: 206202



TRANSDUCERS

Bone Conduction Oscillator and Accessories

Transducer for delivery of acoustic stimuli by bone conduction.

Part #	Item	Qty
BONE-OSC	Silver Navigator Pro Bone Conduction Kit Kit includes: 570-BCB71-1, 541-BONOSC & 570-BCP333	1

NOTE: For use with Bio-logic Navigator Pro AEP and MASTER II systems.



Part #: BONE-OSC

Bio-logic Navigator Pro/ABAer Transducers Insert Earphones (Standard & Broadband)

Part #	Item	Qty
580-SINABR-008	Insert Earphone, 8' cable and Halo Ear Muffin adapter tubes (for ABAer)	1
580-SINSER-012	Insert Earphone, 12' cable	1



Part #: 580-SINABR-008




Part #: EP-HEADPHONE



Part #: 580-AUDISO



Part #: 541-NAVC04



Part #: 101803 - 101898



Part #: 541-TSTCBL



Part #: 001077



Part #: 541-ABRC10-008



Part #: 541-NAVC10

Headphones

Part #	Item	Qty
EP-HEADPHONE	TDH39 Head Phone	1

NOTE: For use with Bio-logic Navigator Pro AEP and MASTER II systems.

Sound Field Accessories

Allows connection of amplified sound field speakers to Navigator Pro hardware.

Part #	Item	Qty
580-AUDISO	Navigator Pro Audio Isolator	1
580-CABFF1	Navigator Pro to Audio Isolator Cable	1
580-AUDISO-NAV	Navigator Pro Audio Isolator Kit	1
	Kit includes: 580-AUDISO & 580-CABFF1	

Bio-logic Navigator Pro — EABR Kit and Parts

Used for performing electrical ABR's.

Part #	Item	Qty
EABR-KIT	EABR Kit Includes:	1
	541-NAVC04, 580-RFFLTR, 3x electrode jumpers	
541-NAVC04	Navigator Pro Trigger In/Out Cable	1
580-RFFLTR	RF Filter Accessory for Navigator Pro	1

Bio-logic Navigator Pro — Electroneuronography Replacement Components

Components connect the Navigator Pro to the Digitimer for electrical stimulation.

Part #	Item	Qty
580-STRB04	Interface Cable, Digitimer BNC to Nav Pro Trigger Cable	1
531-ADAP13	3.5 mm female - 6.35 mm male adapter/ENoG Trigger Cable	1
101803	Reusable Bar electrode - 1 m lead	1
101898	T-Handle for Bar Electrode	1

Bio-logic Navigator Pro — ABaer Hardware Check

Used to perform hardware self diagnostic test for Navigator Pro AEP, MASTER II & ABaer systems.

Part #	Item	Qty
541-TSTCBL	Loop Test Cable	1
001077	Loop Test Cable to Patient Cable Adapter	1

PATIENT CABLES

Cable interface between electrode lead wires from patient to evoked response system.

Part #	Item	Qty
541-ABRC10-008	ABaer Patient Cable, 8'	1
541-NAVC10	Navigator Pro Patient Cable, 12'	1
541-NAVC10-025	Navigator Pro Patient Cable, 25'	1

MISCELLANEOUS ACCESSORIES & SUPPLIES

PC Communication Cables

Part #	Item	Qty
541-USB001	USB A/B Cable, 6' Cable (For NavPro, Scout, AuDX)	1
541-NAVC01	Serial Cable (For ABaer to Panel PC)	1



Part #: 541-NAVC01

Labels for Seiko Label Printer

Part #	Item	Qty
501881	Labels, Seiko Label Printer	1 roll/box
501884	Labels, Seiko Label Printer	4 rolls/box



Part #: 501881

Seiko Label Printer for AuDX, AuDX Pro, & ABaer Systems

Part #	Item	Qty
580-AX-PRINTER	Printer, US Power Supply & Seiko-to-AuDX Serial Cable	1
580-AX-PRINTERE	Printer, EU Power Supply & Seiko-to-AuDX Serial Cable	1
580-AX-PRINTERU	Printer, UK Power Supply & Seiko-to-AuDX Serial Cable	1
541-SEKBLE	Serial Cable, Seiko Printer to AuDX	1
013094	Printer, US Power Supply	1



Part #: 580-AX-PRINTER

Miscellaneous/Other

Part #	Item	Qty
610-CASE08	AuDX/AuDX Pro Standard Carry Case	1
580-AX-PRNCASE	AuDX/AuDX Pro Large Carrying Case (accommodates Printer)	1
CART-AUD-NOTEBOOK	Portable Cart for Notebook Computer - NavPro/ABaer/Scout Sport	1
DX-PRINTER	HP Deskjet® PC Printer for ABaer, Navigator Pro, Scout	1
040947	Barcode Reader (for ABaer)	1
577-NBSCTY	Notebook Computer Lock	1
520-SWBXF	6-outlet Isolation Transformer	1



Part #: 040947

Power Supplies & Cables — Bio-logic Hearing Products

Part #	Item	Qty
002677	Power Supply, 6V DC (For Navigator Pro, AuDX, AuDX Pro, ABaer)	1
003225	Power Cord, 79", US	1
540-LINECD-012	Power Cord, 12", US	1
541-CORD18	Power Cord, 18", US	1



natus[®]

hearing

Natus Medical Incorporated
1501 Industrial Road
San Carlos, CA 94070 USA
1-800-303-0306 • +1-847-949-8615
www.natus.com

© 2015 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.

P/N 015149B