

System requirements

1.5 GHz processor or higher (2 GHz recommended)

2 GB RAM (4 GB recommended)

8 GB free disk space for installation of the OTosuite software.

Microsoft® Windows® 7 SP1 (32 or 64 bit) or Microsoft Windows 8.1 (64 bit) or Microsoft Windows 10 (32 or 64 bit)

All Microsoft Windows updates from Microsoft Update (or .NET 4.6 installed)

USB port for connecting accessories, v. 2.0 or higher

DVD drive (only for a minority of features)

32 bit color display, 1024x768 screen resolution

For Noah mode operation: NOAH 3/NOAH 3 for ENTs (version 3.7 or higher), or Noah 4 (version 4.1 or higher), or OTObase®.

If you are going to use speech material supplied by GN Otometrics, check the Speech Material Installation CD to make sure that your PC has the additional disk space required.