



# Universal programming

## Efficient fitting



### Six ways Hi-Pro 2 boosts fitting efficiency

- Universal programming interface for all hearing instrument brands
- USB connectivity – always ready
- Faster than Hi-Pro predecessors
- No training required
- Sleek, space-saving design
- Backward compatible – works with any hearing instrument that has been programmed by a Hi-Pro device, regardless of version

Hi-Pro 2



**otometrics**  
a division of natus

# Hi-Pro 2 – your connection to greater fitting efficiency



## Streamline your fitting processes

Hi-Pro 2 gives you an easy way to streamline your fitting processes. Developed in collaboration with hearing aid manufacturers, Hi-Pro 2 is the latest in the Hi-Pro series – long recognized as the industry standard for hearing instrument programming. Hi-Pro 2 has a universal hardware interface to program all brands and types of hearing aids – making hearing instrument programming fast, easy and reliable.

## Save time

Hi-Pro 2 saves you time in the fitting process. Programming is faster than its predecessor. The high-speed mode enables fitting modules to communicate with the hearing instrument up to five times faster than the previous Hi-Pro USB version.

### Technical specifications:

#### PC Interface:

Communication: USB 2.0 full speed (USB 1.1 compatible)  
USB Connector Type: "Type B" connector (on the Hi-Pro 2 unit)

#### Power Supply:

The Hi-Pro 2 unit is powered from the PC USB port

#### Operating systems:

Windows XP Pro SP3, Windows 7 (32 and 64 bit),  
Windows 8 (32 and 64 bit), Windows 10

And unlike battery-driven programmers, the Hi-Pro 2 has plug-and-play USB connectivity, so it's ready – and working – when you need it. No training required.

## Boost efficiency

One programming device covers your programming needs. Hi-Pro 2 programs all major hearing instruments. Sleek and discrete, Hi-Pro 2 comes as a space-saving stand-alone device or as an integrated part of the audiometer, Aurical Aud. Hi-Pro 2 is also "backward" compatible, so it works with any hearing instrument that has been previously programmed by a Hi-Pro device, regardless of version.

**Choose Hi-Pro 2 for  
universal programming  
and efficient fitting!**



[Learn more at otometrics.com/hipro](http://otometrics.com/hipro)

**natus.**

**Natus Medical Denmark ApS.** +45 45 75 55 55. [otoinfo@natus.com](mailto:otoinfo@natus.com)  
**Natus Medical Incorporated.** +1 855 283 7978. [otoorder.us@natus.com](mailto:otoorder.us@natus.com)  
[www.otometrics.com](http://www.otometrics.com)

  
**otometrics**  
a division of natus.