

Otometrics® Genie GN200

Space-saving,
high performance
sound booth



Introducing the GN200 single audiometric booth
with optional integrated test equipment



otometrics
a division of natus

The perfect sound booth where space is limited

Designed for single occupancy audiometric testing, the new GN200 offers a high quality-testing environment with exceptional ambient noise control. It is perfect for hearing aid dispensing and testing facilities where space is limited.



The GN200 also comes with an optional Otometrics® Madsen A450 wall mounted audiometer providing you with a ready-to-use test facility with minimal installation required.

You can also combine Video Otoscopy and Tympanometry to deliver a unique integrated audiology solution for you and your client.

All-in-one, integrated sound room

The GN200 is part of the high performance Genie range of sound rooms combined with market-leading Otometrics equipment providing an all in one solution, which is ready to use with minimal set up time.

Easy to install and move

The system is mobile, easy to install and move to new premises if required. It can be shipped flat or fully assembled.

Make it your own

You can also choose your own external branding if required with our bespoke customisation service allowing the exterior of the booth to fit into your clinical or retail environment.

Equipment Options



Otometrics® Madsen A450
Audiometer



Madsen® Zodiac
PC-based Tympanometer



Aurical® Otocam 300
Video Otoscope



Flexible terms

If a single occupancy booth is required as a temporary measure, you can also choose to rent, lease or purchase the system outright.

Inner dimensions

60.9 cm W x 81.3 cm L x 185.4 cm H

Outer dimensions

71.1 cm W x 96.5 cm L x 200.7 cm H

Weight

262 kg

Technical Data

Frequency (Hz)	dB (reduction)
125	23.3
250	28.7
500	34.8
1000	37
2000	39.6
4000	46.3
8000	x

Ref: Mecanum Tests Report MNS20170816-01

Features

- Integrated adjustable high comfort padded seat easy to clean and maintain
- Integrated jack panel with USB connectivity
- Rood mounted low profile HALO LED lighting (CE approved)
- Integral ventilation system
- High performance acoustic glazing
- Magnetically sealed door with recessed handle. Door may be rotated to suit left or right opening
- Audiometer shelf with brackets or wall mountable audiometry options
- High quality paint finish with a choice of colours
- Optional external wrapping service to allow for unique branding to the individual clinic or retail environment
- Easy to install and move. Installation service available and may be delivered fully assembled or flat packed

Learn more

Call our office on **0870 9000 675** or come and see our showroom at the following address:

Natus Nicolet UK Ltd.

Unit 4, Cambridge Terrace, St James Road, Brackley
NN13 7XY Northamptonshire

You can also email us at otoinfo.uk@natus.com for more information.

About Otometrics

Otometrics, a division of Natus, is the audiology industry leader providing instrumentation and software solutions to hearing and balance care professionals worldwide. Otometrics help hearing and balance care professionals improve the quality of life for clients and patients by delivering expert knowledge, reliable solutions and trusted partnership. Otometrics is a strategic business unit of Natus Medical Incorporated (North America NASDAQ: BABY). For more information, please visit www.otometrics.com or www.natus.com.

natus

Natus Medical Denmark ApS. +45 45 75 55 55. otoinfo@natus.com
Natus Nicolet UK Ltd. +44 870 9000 675. otoinfo.uk@natus.com
www.otometrics.co.uk



otometrics
a division of natus