



## Ensure consistency and quality in every fitting

### Required Measurements

Smooth fitting processes in a consistent quality are key to ensuring your clinic or practice runs efficiently and with the right level of quality and service. Now you can ensure that everyone who fits clients in your practice performs the audiometry and verification measures you define as a standard part of your protocol – without disrupting their personally preferred workflow. Required Measurements in Aurical helps you to streamline fitting processes across your entire practice even more.



Required Measurements is a dedicated feature in Aurical available in the audiometry and PMM modules in Otosuite®. This enables you to identify the measurements that must be included during the visit and allows saving of the data only after these measurements are completed. If a required measurement has not been performed, a dialogue box notifies the user and displays the tests that need to be completed. If the measurement cannot be completed, the user can enter an explanation in a text box, and the results can be saved.



# Consistency and quality with every fitting

## Required Measurements offers:

- **Easy configuration** You or your clinic administrator define the measurements you want to be performed in the Required Measurements Editor and deploy a simple file to each Otosuite user.
- **Quick user notification and a clear overview of the data that is missing** Users receive a reminder on their screen if Audiometry or Probe Microphone Measurements are missing at the end of their session.
- **Simple comment field for user feedback and communication** Users can provide feedback regarding missing data via an easy-to-use comment field.

Aurical, Otosuite and Required Measurements give you yet another way to streamline your fitting processes while ensuring quality and consistency with every fit.

Learn more about our complete fitting solution. Contact your local Natus representative today or visit [hearing-balance.natus.com/aurical](http://hearing-balance.natus.com/aurical)



A reminder appears on the screen if users are missing any data at the end of their session.

A complete, modular fitting solution, Aurical encompasses audiometry, PMM, HIT, video otoscopy and a HI programming interface in one system. Hearing care professionals can use the modules independently or combine them as needed. This latest release demonstrates once again how Aurical continues to innovate – enabling hearing care professionals to realize greater efficiency gains while boosting the quality of counseling and care.



Healthcare solutions with one thing in mind. You.

©2020 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates. 7-26-1670-EN Rev02

**natus**®

Natus Medical Incorporated

[natus.com](http://natus.com)