




## Single wall audiometric sound rooms

Model number	Dimensions	
	Interior	Exterior
<b>GN-11</b>	2'9" x 3'8" x 6'6"H (83.8 cm x 111.8 cm x 198.1 cm)	3'5" x 4'4" x 7'6" (104.1 cm x 132.1 cm x 228.6 cm)
<b>GN-12</b>	3'8" x 4'5" x 6'6"H (111.8 cm x 139.7 cm x 198.1 cm)	4'4" x 5'3" x 7'6" (132.1 cm x 160.0 cm x 228.6 cm)
<b>GN-12A</b>	3'8" x 7'4" x 6'6" (111.8 cm x 223.5 cm x 198.1 cm)	4'4" x 8'0" x 7'6" (132.1 cm x 243.8 cm x 228.6 cm)
<b>GN-13</b>	4'7" x 5'6" x 6'6" (139.7 cm x 167.6 cm x 198.1 cm)	5'3" x 6'2" x 7'6" (160.0 cm x 188.0 cm x 228.6 cm)
<b>GN-13A</b>	4'7" x 7'4" x 6'6" (139.7 cm x 223.5 cm x 198.1 cm)	5'3" x 8'0" x 7'6" (160.0 cm x 243.8 cm x 228.6 cm)
<b>GN-14</b>	5'6" x 6'5" x 6'6" (167.6 cm x 195.6 cm x 198.1 cm)	6'2" x 7'1" x 7'6" (188.0 cm x 215.9 cm x 228.6 cm)
<b>GN-15</b>	6'5" x 7'4" x 6'6" (195.6 cm x 223.5 cm x 198.1 cm)	7'1" x 8'0" x 7'6" (215.9 cm x 243.8 cm x 228.6 cm)
<b>GN-16</b>	7'4" x 8'3" x 6'6" (223.5 cm x 251.5 cm x 198.1 cm)	8'0" x 8'11" x 7'6" (243.8 cm x 271.8 cm x 228.6 cm)
<b>GN-17</b>	8'3" x 9'2" x 6'6" (251.5 cm x 279.4 cm x 198.1 cm)	8'11" x 9'10" x 7'6" (271.8 cm x 299.7 cm x 228.6 cm)
<b>Corner designs</b>		
<b>GN-46</b>	5'10" x 5'10" x 6'6" (177.8 cm x 177.8 cm x 198.1 cm)	6'6" x 6'6" x 7'6" (198.1 cm x 198.1 cm x 228.6 cm)
<b>GN-66</b>	6'6" x 6'6" x 6'6" (198.1 cm x 198.1 cm x 198.1 cm)	7'2" x 7'2" x 7'6" (218.4 cm x 218.4 cm x 228.6 cm)
<b>GN-86</b>	6'6" x 8'6" x 6'6" (198.1 cm x 259.1 cm x 198.1 cm)	7'2" x 9'2" x 7'6" (218.4 cm x 279.4 cm x 228.6 cm)




## Double wall audiometric sound rooms

Model number	Dimensions	
	Interior	Exterior
<b>GN-24</b>	5'6" x 6'5" x 6'6" (167.6 cm x 195.6 cm x 198.1 cm)	7'6" x 8'5" x 8'2" (228.6 cm x 256.5 cm x 248.9 cm)
<b>GN-25</b>	6'5" x 7'4" x 6'6" (195.6 cm x 223.5 cm x 198.1 cm)	8'5" x 9'4" x 8'2" (256.5 cm x 284.5 cm x 248.9 cm)
<b>GN-26</b>	7'4" x 8'3" x 6'6" (223.5 cm x 251.5 cm x 198.1 cm)	9'4" x 10'3" x 8'2" (284.5 cm x 312.4 cm x 248.9 cm)
<b>GN-27</b>	8'3" x 9'2" x 6'6" (251.5 cm x 279.4 cm x 198.1 cm)	10'3" x 11'2" x 8'2" (312.4 cm x 340.4 cm x 248.9 cm)



## Single wall control/single wall exam audiometric suites

Model number	Dimensions		
	Control Room Interior	Exam Room Interior	Exterior
<b>GN-34</b>	5'6" x 4'7" x 7'6" (167.6 cm x 139.7 cm x 198.1 cm)	5'6" x 6'5" x 7'11" (167.6 cm x 195.6 cm x 241.3 cm)	6'2" x 12'8" x 7'6" (188.0 cm x 386.1 cm x 228.6 cm)
<b>GN-45</b>	6'5" x 5'6" x 6'6" (195.6 cm x 167.6 cm x 198.1 cm)	6'5" x 7'4" x 7'11" (195.6 cm x 223.5 cm x 241.3 cm)	7'1" x 12'8" x 7'6" (215.9 cm x 442.9 cm x 228.6 cm)
<b>GN-56</b>	7'4" x 6'5" x 6'6" (223.5 cm x 195.6 cm x 198.1 cm)	7'4" x 8'3" x 7'11" (223.5 cm x 251.5 cm x 241.3 cm)	8'0" x 14'6" x 7'6" (215.9 cm x 497.8 cm x 228.6 cm)
<b>GN-67</b>	8'3" x 7'4" x 6'6" (251.5 cm x 223.5 cm x 198.1 cm)	8'3" x 9'2" x 7'11" (251.5 cm x 279.4 cm x 241.3 cm)	8'11" x 18'2" x 7'6" (271.8 cm x 553.7 cm x 228.6 cm)



## Single wall control/double wall exam audiometric suites

Model number	Dimensions		
	Control Room Interior	Exam Room Interior	Exterior
<b>GN-121</b>	7'4" x 5'6" x 7'9" (223.5 cm x 167.6 cm x 236.2 cm)	6'0" x 6'8" x 6'6" (182.9 cm x 203.2 cm x 198.1 cm)	8'0" x 14'6" x 8'2" (243.8 cm x 442.0 cm x 248.9 cm)
<b>GN-122</b>	6'5" x 8'3" x 7'9" (195.6 cm x 251.5 cm x 236.2 cm)	6'11" x 7'7" x 6'6" (210.8 cm x 231.1 cm x 198.1 cm)	8'11" x 16'4" x 8'2" (271.8 cm x 497.0 cm x 248.9 cm)
<b>GN-123</b>	9'2" x 7'4" x 7'9" (279.4 cm x 223.5 cm x 236.2 cm)	7'10" x 8'6" x 6'6" (238.8 cm x 259.1 cm x 198.1 cm)	9'10" x 18'2" x 8'2" (299.7 cm x 553.7 cm x 248.9 cm)
<b>GN-124</b>	10'1" x 8'3" x 7'9" (307.3 cm x 251.5 cm x 236.2 cm)	8'9" x 9'5" x 6'6" (266.7 cm x 287.0 cm x 198.1 cm)	10'9" x 20'0" x 8'2" (327.7 cm x 609.6 cm x 248.9 cm)

