



Feel the difference

Quick and efficient immittance screening




MADSEN® Zodiac is an immittance testing solution from Otometrics. Built from the ground up with hearing professionals like you, MADSEN Zodiac delivers the control, confidence and efficiency you need for successful immittance testing.

MADSEN Zodiac is designed with focus on the probes and the way you work. It is reliable, responsive and easy to use. You stay focused on your patients, boost testing efficiency and increase testing success. MADSEN Zodiac is suitable for use with both children and adults.

MADSEN® Zodiac

 **otometrics**




Test more patients every day

The Quick Check version of MADSEN® Zodiac makes tympanometry and reflex screening more efficient. You can test more patients with greater control every day – or easily introduce immittance testing as part of your hearing care protocol for a more efficient workflow. The design is clean and simple – perfect for screening – with only the buttons and essential tests you need. The flat menu structure with lighted, one touch/one function buttons provides easy navigation, allowing you to focus on the patient. The PC-based version runs on the intuitive OTOsuite® software.

Easier way to a reliable seal

The Screener probe has a unique ergonomic design for easier, faster and more stable placement. It feels natural in your hands and you can hold it in the way that works best for you to make it easier to place the probe and start testing.

Now it is no longer necessary to walk back and forth to the device. When you reach the seal, testing starts automatically. You can change ears whenever you want with the remote control button on the probe. And long, sturdy cables allow you to move around freely and get closer to your patient.



No-hassle print and easy-to-change paper roll. The built-in printer prints test results and reports quickly. OTOsuite software gives you the choice of more printing options.

MADSEN Zodiac streamlines immittance screening. Whether you choose the PC-based or stand-alone version, MADSEN Zodiac is simple to use. Just walk up and start testing – it really is that easy and intuitive.

Contact your local Otometrics representative for a demonstration. Or, visit www.otometrics.com/zodiac for more information.