



Patient database system

The OTObase is a patient management solution providing a user-friendly customer interface for efficient data management. The system works in two configurations; with a patient management system (PMS) and as a stand-alone system. Data storage can happen in the database that is an integrated part of the product.

The graphical user interface provides a fast overview of both the measurement data and the stored Muster 15 reports for the selected patient. The data and the report storage are presented as sessions with clear indication of time and date. Furthermore, the data is easily accessed by clicking the data entry. Hereby either the measurement module or the muster 15 report comes up.

Patients are added/deleted or opened for editing with a single click. Instant search provides a fast way to selecting the patient. The graphical user interface is clean and structured to give a fast overview of the data for the individual patient

Intuitive patient management

Password protection

Import of Noah data

SQL database included

Works with own database

Easy access to external measurement modules

List of measurement session

GDT interface for data transfer with Patient Management Systems

Two page Muster 15 report

Tone and speech data from different sessions

List of created reports



Technical Specifications	
Supported OS:	Version 1.3 supports Windows 7/8/10 32/64 bit
NOAH supported:	NOAH 4.0
SQL servers supported:	2005 and 2008
Muster 15 report:	Yes
GDT Interface:	Yes

About Otometrics

Otometrics is the audiology industry leader providing instrumentation and software solutions to hearing and balance care professionals worldwide. For more than 50 years, Otometrics has been helping hearing and balance care professionals succeed in improving the quality of life for their clients and patients around the world by delivering expert knowledge, reliable solutions and services – providing a trusted partnership.

Otometrics develops, manufactures and markets computer-based audiological, otoneurologic and vestibular instrumentation in more than 80 countries under the Otometrics brands, including MADSEN®, AURICAL®, HORTMANN®, ICS® and Oscilla. Otometrics is part of Natus Medical Incorporated (North America NASDAQ: BABY) – a US-based leading provider of healthcare products, devices and services used in the newborn, neurology, sleep, hearing and balance markets. For more information, please visit www.otometrics.com and www.natus.com.

facebook.com/otometrics

twitter.com/otometrics