



## Document Description

The document describes the function of the “Protect User Test” setting.


1. In OTOsuite it is possible to block a user from saving a user test or a startup setting due to the “**Save user test**” and the “**Save options as start up**” icons.




1. To activate the Icons again, go to “**Tools**” and “**Configuration Wizard**”.




2. Click on **“Next”**.




3. Click on **“Next”**.



4. Click on **“Next”**.



5. Remove the check mark from **“Protect User Tests”**



6. Now the **“Save user test”** and the **“Save options as start up settings”** icons are active.

